**Produktgruppe**

**Product Group**

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## DISCLAIMER

Dear Customer,

Thank you for choosing a MAVIG product. Please read this document carefully and follow all safety relevant instructions and notes.

The manual is part of the product and must therefore be kept in a place close to it so that important information and safety instructions for use can be looked up at any time.

Do not pass the product on to third parties without the valid instructions for use.

Please contact us if you have questions or comments about our products.

## Contact information

If you have any questions or comments, you will find contact information for your local MAVIG subsidiary at the end of this manual.

## Disclaimers

The content of this manual is for information purposes only. The products described are subject to change without notice in the course of continuous development.

MAVIG only assumes warranty for this manual and the products within the scope of the manufacturer's liability.

MAVIG is not responsible for unauthorized modifications made to the product or use of the product for unintended purposes. MAVIG cannot be held liable for improper operation and modifications. Since improper modifications may affect the correct operation, safety or reliability of the product, product modifications require written authorization from MAVIG.

## Changes

MAVIG products are constantly being further developed. MAVIG reserves the right to change the manual, packaging, equipment or technical specifications as well as scope of delivery without notice.

## Translations

The original document is in German. In the event of a translation, the German version shall remain the authoritative version.

## Copyright and trademarks

All rights reserved. These instruction is protected by copyright. MAVIG is a registered trademark. All other trademarks mentioned in these instruction are the property of the respective manufacturer.

## INDEX

<b>1.</b>	<b>General information</b>	<b>13</b>
1.1.	Manufacturer	13
1.2.	Operation and Operator	13
1.3.	Symbols, Safety and Warning instructions	13
1.4.	Reporting Requirements according to Regulation (EU) 2017/745	13
<b>2.</b>	<b>Product information</b>	<b>14</b>
2.1.	Product description	14
2.2.	Intended Use	14
2.3.	Contraindications and Exclusions	14
2.4.	Identification	14
2.5.	CE-Marking	14
2.6.	Accessories and Spare Parts	14
<b>3.</b>	<b>Installation and Use</b>	<b>15</b>
3.1.	User Environment	15
3.2.	Installation, Connection methods, Connections	15
3.3.	Cleaning, Disinfection, Sterilization	15
3.3.1.	Cleaning	15
3.3.2.	Disinfection	16
3.4.	Maintenance and Repair	17
3.5.	Calibration and Verification	17
3.6.	Warranty	17
<b>4.</b>	<b>Storage, Transport, Shelf Life</b>	<b>17</b>
<b>5.</b>	<b>Technical Data</b>	<b>18</b>
<b>6.</b>	<b>Decommissioning and Disposal</b>	<b>18</b>

## 1. General information

### 1.1. Manufacturer

MAVIG GmbH  
Stahlgruberring 5  
81829 Munich  
Germany

Phone +49 (0) 89 420 96 - 0  
Fax +49 (0) 89 420 96 - 200

eMail [info@mavig.com](mailto:info@mavig.com)  
Web [www.mavig.com](http://www.mavig.com)

### 1.2. Operation and Operator

Only appropriately trained personnel may operate, clean and disinfect the product (medical and radiological technical knowledge in accordance with country specific regulations).

The instructions for use are valid only after a properly conducted first commissioning.

Product must be thoroughly cleaned and disinfected before first use.

### 1.3. Symbols, Safety and Warning instructions

<b>SAFETY-RELEVANT</b>	Safety relevant instructions are mandatory and must be followed. Non-observance could result in <b>injury</b> or <b>death</b> of one or more persons. You must read and understand these instructions before starting up the systems.
<b>CAUTION</b>	Cautions must be observed, otherwise the device may be <b>damaged</b> or <b>malfunction</b> . You must read and understand these instructions before starting up the system.
	<b>Not sterilizable</b> The product cannot be sterilized with ethylene oxide, irradiation, aseptic processing techniques, steam or dry heat.
	<b>Medical Device</b> Indicates that the product is a medical device.

### 1.4. Reporting Requirements according to Regulation (EU) 2017/745

In accordance with the relevant legal regulations, serious incidents that have occurred in connection with a MAVIG medical device must be reported to MAVIG and to the responsible national authorities in the country in which the user is active. In this sense, incidents are reportable if the use of the medical device directly or indirectly led, might have let or might lead to the death or the temporary or permanent serious deterioration of a patient's, user's or other person's state of health

## 2. Product information

### 2.1. Product description

Please inform yourself about product variants and technical performance in our product data sheets and product catalogs on our website [www.mavig.com](http://www.mavig.com).

### 2.2. Intended Use

Radiation protective shields are used for the stationary protection of medical personnel against ionizing radiation. The product is designed for applications in air-conditioned health care facilities such as clinics, hospitals and doctors' surgeries. The product is used exclusively by medical professionals such as physicians and assistants.

### 2.3. Contraindications and Exclusions

Any other purpose than the one described is not intended. The manufacturer is not liable for any resulting damage.  
Do not install in MRI rooms.

### 2.4. Identification

The product can be clearly identified by using the type plate / product label attached on the product. There you will find at least the following information:  
*name of the product, manufacturer and address, serial number*

### 2.5. CE-Marking

The products meet the requirements of the EU Medical Device Regulation 2017/745 (MDR). The CE conformity of the product is identified by the CE symbol on the label.

### 2.6. Accessories and Spare Parts

Only the original manufacturer's components and spare parts must be used with respect to this product. Failure to comply with this instruction will exempt the manufacturer from any liability or warranty obligations. The operator undertakes to ensure that the device is used only with permitted replacement parts. The use of unauthorized parts entails unknown risks and must be avoided at all times.

For information about available spare parts and accessories, please refer to our data sheets, which are available at our website [www.mavig.com](http://www.mavig.com) or contact us at:

eMail [info@mavig.com](mailto:info@mavig.com)      oder      Telefon      +49 (0) 89 420 96 – 0

#### **SAFETY- RELEVANT**

Use of spare parts from third-parties can impact the safety and result in damage, incorrect function or total failure of the product. Only use original spare parts!

### 3. Installation and Use

#### 3.1. User Environment

During the use of this product, the room temperature and humidity should remain in a range for normal working and living conditions.

#### 3.2. Installation, Connection methods, Connections

The MAVIG product/ system must be set up correctly in accordance with the assembly instruction enclosed to the product/ system. Failure to follow the instruction can cause the system to be unstable and/or mean that its function is reduced.

The precise installation is explained in document POR03Oxx

#### 3.3. Cleaning, Disinfection, Sterilization


##### 3.3.1. Cleaning

For wiping all surfaces:

- Use soft, particle-free and free of residues of other cleaning agents cloths for cleaning.
- Only slightly moisten the cloth. If necessary, use a mild soap solution or commercial washing-up liquid.
- Make sure that no liquid is spilled during cleaning and penetrates the device.
- Then rub dry the surfaces with a cloth.

<b>CAUTION</b>	<p><b>Damage to the product</b></p> <p>In order to avoid damage to the product no aggressive, abrasive, caustic, bleaching, alkaline, acidic or corrosive agents should be used (e.g. solvents, acetates, scouring agents, chlorine etc.).</p>
----------------	--

Only valid for OT54xx sowie OT94xx: Cleaning of ComforTex®HPMF Curtain  
The following markings apply to the X-Ray protection curtain/strips:

Fresh stains can usually be removed from the radiation protective clothing with soapy water of approx. 20°C and a soft cloth. For stubborn stains, we recommend a 10% solution of sodium thiosulphate in water. For fluids containing proteins (such as blood), use cold water. Please make sure that the cloth is damp and not wet.

### 3.3.2. Disinfection

Only use wipe disinfection.

- No liquid may penetrate the product.
- The applied reprocessing process must comply with the applicable legal regulations and guidelines.
- Hygiene procedures and conditions are the responsibility of the operator and must be determined by him.
- Ensure that the materials used are compatible (e.g. stainless steel, aluminum, plastic, lead acrylic, lead glass, PVC, natural, silicone or synthetic rubber, chromed, painted or coated surfaces).
- In case of doubt, test the product suitability in an inconspicuous place.

<b>CAUTION</b>	<b>Surface Damage</b> The combined use of different disinfectants and/or disinfectants and cleaning agents may impair the product performance or damage surfaces of the product.
----------------	---

Only valid for OT54xx sowie OT94xx : Disinfection ComforTex®HPMF Curtain

Always perform cleaning prior disinfection.

<b>CAUTION</b>	<b>Damage of Product</b> Do not use alcohol-based disinfectants, pure alcohol (e.g. ethanol and isopropanol), or oxidizing as well as active chlorine or salt-containing disinfectants as they attack the surface of the outer material more strongly than disinfectants based on other chemicals.
----------------	---

Avoid applying too much pressure while disinfecting. To maintain the properties of our radiation protection products, we recommend the use of the following disinfectants. Please observe the instructions below:

- The tests carried out at our company do not correspond to any certified or standardized test procedures, so that we cannot issue a release, but only a recommendation.
- Since we have no influence on the production as well as the raw materials of the disinfectant, we have to decline any guarantee on the resistance of our protective products to the disinfectant.
- The effectiveness of the surface disinfectant has not been considered. The instructions for use of the disinfectant must be adhered to.
- For the reasons mentioned above MAVIG excludes any warranty or liability for subsequent damages.

Disinfectant	Manufacturer	Effective base
Microbac forte	BODE Chemie GmbH	Quaternary compounds, alkylamine derivative
Kohrsolin FF	BODE Chemie GmbH	Aldehyde, quaternary compounds
Antifect FF	Schülke & Mayr GmbH	Quaternary compound, aldehyde
Microzid Sensitive Wipes	Schülke & Mayr GmbH	Quaternary compounds
Microbac Tissues Solution	BODE Chemie	Quaternary compounds

In general, when choosing a disinfectant, make sure it contains as few alcohols as possible. Ideally, the first 5 ingredients should not be alcohols or these should be listed with amounts of „< 5%“.

### 3.4. Maintenance and Repair

MAVIG GmbH distributes its products through authorized dealers only. Those dealers possess the necessary knowledge for servicing, maintenance and repair of the product. Please approach your authorized dealer first in case you need any of the above mentioned services.

Instruction for frequency and coverage of maintenance are to be found in document POR02Gxx.

### 3.5 Calibration and Verification

<b>CAUTION</b>	Check for visible damage prior to usage of product. Do not use a damaged product!
----------------	---

### 3.6. Warranty

Our general terms of business apply. These should be available to the operator latest after completion of the sales contract.

MAVIG guarantees the safety and reliability of the product/ system solely on the following basis:

- The product/system is used solely for its intended purpose and its use complies with the information provided in these instructions for use.
- Only original authorised replacement parts or accessories as specified by MAVIG are used. The use of unauthorised parts holds unknown risks and must be avoided at all times.
- No unauthorised structural or other changes are to be made to the product/system. Unauthorised changes or modifications to the product are prohibited for safety reasons.
- The product/system is regularly inspected and serviced at the specified times.

The manufacturer's warranty period for workmanship and materials is 24 months.

## 4. Storage, Transport, Shelf Life

Stored only in closed or covered spaces.

Transport and storage conditions:

- Temperature: from -20°C to +70°C
- Relative humidity: max. 95%
- Air pressure: from 700 hPa to 1060 hPa

In service:

- Temperature: from +10°C to +70°C
- Relative humidity: from 20% to max. 75%
- Air pressure: from 700 hPa to 1060 hPa

## 5. Technical Data

Please inform yourself about the technical performance of the product in our product data sheets and product catalogs, which are published on our website [www.mavig.com](http://www.mavig.com) or contact us at:

eMail	<a href="mailto:info@mavig.com">info@mavig.com</a>	Phone	+49 (0) 89 420 96 - 0
-------	--	-------	-----------------------

## 6. Decommissioning and Disposal

The device must be taken out of operation if it is damaged or has not passed the maintenance test.

The device must be disposed of in accordance with your country-specific regulations.

Please feel free to contact us if you have any questions:

eMail	<a href="mailto:info@mavig.com">info@mavig.com</a>	Phone	+49 (0) 89 420 96 - 0
-------	--	-------	-----------------------

### MAVIG GmbH Headquarters

PO Box 82 03 62  
81803 Munich  
Germany

Stahlgruberring 5  
81829 Munich  
Germany

Phone +49 (0) 89 420 96 0  
Fax +49 (0) 89 420 96 200  
e-Mail [info@mavig.com](mailto:info@mavig.com)

### Nordic & Baltic Countries MAVIG Nordic

Stockholm  
Sweden

Phone +46 (0) 722 25 25 68  
e-Mail [larsson@mavig.com](mailto:larsson@mavig.com)

### Benelux, GB, Ireland MAVIG B.V.

Mercuriusweg 86  
2516 AW Den Haag  
Netherlands

Phone +31 (0) 70 33 11 688  
Mobile +31 (0) 61 595 43 48  
e-Mail [simmonds@mavig.nl](mailto:simmonds@mavig.nl)

### France MAVIG France SARL

66, Ave. des Champs Elysées  
F-75008 Paris  
France

Phone +33 (0)1 30 59 46 23  
Fax +33 (0)1 30 59 46 23  
e-Mail [info@mavig.fr](mailto:info@mavig.fr)

### Southeast Europe MAVIG Southeast Europe

Ljubljana  
Slovenia

Phone +386 (40) 6 33 900  
e-Mail [dejak@mavig.com](mailto:dejak@mavig.com)

### China MAVIG Healthcare Beijing Co., Ltd.

Room #1010, Floor 9  
#609 Wang Jing Yuan  
Chaoyang District  
Beijing 100102, P.R.China

Phone +86 - 10 - 64 78 19 21  
Fax +86 - 10 - 64 78 19 21  
Mobile +86 - 13 70 122 85 86  
e-Mail [zhang@mavig.com](mailto:zhang@mavig.com)

### REPRESENTATION

#### USA & Canada

Ti-Ba Enterprises, Inc.

25 Hytec Circle  
Rochester, NY 14606  
USA

Phone +1 (585) 247 1212  
Fax +1 (585) 247 1395  
e-Mail [mavigusteam@mavig.com](mailto:mavigusteam@mavig.com)

#### Middle East & Northern Africa Mena Medical Development

Ashrafieh – Adlieh Square – Alfaras St.  
- Alboustany Building, 5th floor  
Beirut  
Lebanon

Phone +961 14 23 499  
Fax +961 14 26 499  
e-Mail [mavig@mena-md.com](mailto:mavig@mena-md.com)

#### Russian Federation

AO Sante Medical Systems

Novodmitrovskaya 2, bld.1  
127015 Moscow  
Russia

Phone +7 499 551 55 73/75  
e-Mail [info@sante.ru](mailto:info@sante.ru)

**[www.mavig.com](http://www.mavig.com)**