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For MAVIG Systems

List of contents

1. Safety instructions	page 4
2. Brief description of the light YLED-1F	page 8
3. Operating the light YLED 1F	page 9
3.1 ON/OFF switch	page 9
3.2 Light intensity control	page 9
3.3 Endo-light	page 9
3.4 Focusing	page 10
3.5 Positioning	page 10
3.6 Addition to positioning the lamp	page 11
3.7 Danger of collision while positioning the lights	page 11
4. Cleaning	page 12
4.1 Sterilisable handle	page 12
4.2 Lamp head, lens and protective disk	page 14
5. Maintenance	page 15
5.1 Regular maintenance work	page 15
6. Data	page 16
6.1 Technical data	page 16
6.2 Electrical data	page 17
6.3 Information regarding the electrical installation	page 17
6.4 Weights	page 17
6.5 Environmental conditions	page 18
6.6 General remarks	page 19
7. CE-mark	page 19
8. Manufacturer	page 19
9. Disposal	page 19
10. Electromagnetic compatibility	page 20
11. Overview of the light functions	page 26

Dear customer!

Congratulations for acquiring our OT-light **YLED-1F**.

The new OT-light generation with LED technology supports your professionalism by innovative technology and design.

The advantages of the LED technology: a life-span of 60.000 hours and an almost non-existent heat development in the surgeon's head area and in the wound field.

The advantages already provided by Dr. Mach's light technology with halogen and gas discharge lamps have been maintained: natural colour reproduction, exact illumination of the wound field, and easy positioning of the light head.

For detailed information for Installation, please use the MAVIG installation manual for Portegra2 360°, or TWIN or Ondal Acrobat Swing.

The information in this manual has been based on information from the **manufacturer Dr. Mach GmbH & Co. KG, Germany**
MAVIG does not assume any liability for this information.

For questions please contact:

MAVIG GmbH, Stahlgruberring 5, 81829 Munich
Tel.: +49 (0) 89 42096-0, Fax +49 (0)89 42096-200
Internet: www.mavig.com, E-Mail: info@mavig.com

1. Safety instructions

Pay attention to the instructions for use when handling the lamp.

WARNING:

This device has not been designed for use in potentially explosive areas.
According to the Medical Device Regulation the light is classified as class I.

Transport damage:

When you receive the goods, please check them for transport damages.
Transport damages must be reported to the supplier immediately.

Store the light in its package for at least 24 hours in the respective room before mounting, in order to equal temperature differences.

Please read the instructions for use carefully to make the most of your lighting system and to avoid any damages to the device.

The light may only be repaired and special assembly work may only be carried out by ourselves or a company that has been expressly authorized by us.

The manufacturer can only be made responsible for the safety of the light if repairs and alterations are carried out by the manufacturer himself or a company that guarantees to observe the safety regulations.

No modification of the lamp is allowed!

The manufacturer cannot be made liable for personal or material damages if the light is operated inexpediently or incorrectly or used for purposes other than those for which it is intended.

The light is to be dismantled from the spring arm in reverse order to its assembly. This may only be carried out after the spring arm has been adjusted in height at horizontal position since the arm is under spring tension and can bounce up.

Make sure that the light is in perfect working order before every use.

Attention!

A main control switch must be installed for turning the system power-off.

ATTENTION!

During the mounting of the light the entire system (incl. the ceiling attachment) must be disconnected from mains!
A later dismantling of the light from the spring arms or dismantling the sliding contacts inside the arms is to be done **ONLY AFTER DISCONNECTING THE ENTIRE SYSTEM FROM MAINS.**
Otherwise the electronic board will be damaged!

Symbols and notes used in this user manual:


This symbol means possible hazard sources. Please observe also the safety remarks and the hazard specifications mentioned in the mounting instructions and user manuals from MAVIG company.

This symbol means possible hazard caused by electric current. Please observe also the safety remarks and the hazard specifications mentioned in the mounting instructions and user manuals from MAVIG company.

This symbol refers to important mounting indications, useful information and operation hints.

Temperature range for transport and storage

Indication for disposal

CE- conformity mark

This way up

Keep dry

Fragile, Handle with care

Symbols and notes used on the device:


This symbol indicates to consult the user manual.

This symbol indicates that the user manual is digital available in further languages.

Indication on China RoHS / Pollution control Logo China

Temperature range for transport and storage

Indication for disposal

Serial number of the product

Article number of the product

Address of manufacturer / distributor of the product

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2. Brief description of the light YLED-1F

YLED-1F intended use:

The YLED-1F lighting system is designed for illuminating an examination area at the hospital and doctor's practice.

YLED-1F indications for use:

The YLED-1F is intended to illuminate the surgical field and the patient's body with a high intensity, shadow-free and "cold" light.

Essential Performance:

The examination light YLED-1F is intended to provide the illumination.

General product description

- The YLED-1F lighting system is an examination light according to EN 60601-2-41, which is not fail-safe when used as a single light.
- The YLED-1F lighting system is designed to support therapy and diagnosis.
- The light is used in medical rooms (groups 0, 1 and 2 according to DIN VDE 0100-710).
- This light system can be added to the ceiling mounted suspension system supporting the horizontal arms and spring arms, as well as a wall light or mobile light.
- The maintenance of the light must be done every two years.
- The electrical connection for the ceiling and wall lights is done by a fixed connection.

Accessories:


The accessory for the YLED-1F lighting systems is as follows:

- Sterilizable handle

Spare Parts List:

- Sterilizable Handle	Article Code:	E-7110081
- Clear protective Cover	Article Code:	E-7110981
- Control Panel	Article Code:	E-7110982
- Upper Lamp Housing	Article Code:	E-7110983

3. Operating the light YLED-1F

3.1 ON/OFF switch

The push button **1** on the control panel turns the light YLED-1F ON and OFF.

3.2 Light intensity control

The light YLED-1F offers the facility of light intensity control.


The adjustment range of the light intensity is from 50 % to 100 %.

The light intensity can be adjusted according to the requirements of the surgeon / physician.

The light intensity can be increased by pressing push button **2**.

The light intensity can be decreased by pressing push button **3**. A further press of button **3** activates the endo-light function of the lamp.

The set light intensity is shown by LEDs **4**. One illuminated LED indicates the lowest possible intensity setting; five illuminated LEDs indicate the highest possible intensity setting.

3.3 Endo – Light

In case of endoscopic intervention the light of the examination light can be dimmed.

The Endo mode is activated by repeated pressing push button **3**, until only LED **4** is on. When pressing the button again, LED **4** goes off and LED **5** turns on automatically.

The Endo mode is now active.

This mode can be deactivated by pressing the button **2**.

3.4 Focusing

The lamp-models YLED-1F have a focusing function. That means, you can either enlarge the diameter of the light field or bundle the light to a smaller area, depending on the circumstances.


To activate the function of focusing turn the handle **6** (see figure).

3.5 Positioning

Use the handle **6** to position the lamp.

Use the handle for positioning the light during the operation.

(Example Picture)


3.6 Addition to positioning the lamp

When the spring arm and extension arm are in a 180° position, the lamp head cannot be moved in the direction of the suspension axis.

A repeated force of > 25 N (according to DIN EN 60601-2-41) at this angle can lead to damage.

This applies for all ceiling and wall mounted lamps (example in picture on the left).

In order to be able to move the lamp in the direction of its suspension axis, the angle between the extension arm and spring arm must be < 180°.


3.7 Danger of collision while positioning the lights

During positioning, eventual collisions between the lights, spring arms and other devices must be avoided.

Cover parts can get loose and fall down.

4. Cleaning

4.1 Sterilizable handle

The light is equipped with the sterilisable handle **6**. The handle sleeve is removable and sterilisable. Before using the first time and before every use the handle sleeve must be cleaned, disinfected and sterilised.

The handle sleeve must be removed for sterilisation:

- To remove press the lock **V** and pull off the sterilizable handle sleeve **6** while keeping the lock pressed.
- To attach, push on and slightly twist the handle until the lock **V** engages securely.

Handles will become unsterile during an operation. Therefore always keep additional handles available for exchange.

Cleaning / disinfection and sterilisation

Basics

Efficient cleaning / disinfection is an essential requirement for effective sterilisation of the handle.

Within the scope of responsibility for the sterility of the products it should be noted that only sufficiently validated equipment and product specific processes are used for cleaning / disinfection and that the validated parameters are complied with in every cycle.

In addition, the hospital / clinic hygiene regulations must be observed.

Cleaning / disinfection

Cleaning and disinfection must be carried out immediately after use.

A mechanised process (disinfector) should be used for cleaning / disinfection. The efficiency of the process used must be recognised and validated in principle (e.g. listed under disinfectants and disinfection procedures tested and recognised by Robert-Koch-Institute / DGHM).

When using other procedures (e.g. a manual procedure), proof and process efficiency in principle must be provided within the scope of validation.

Proof in principle of the suitability of the handles for efficient cleaning / disinfection was provided using a cyclic cleaning system (Netsch-Bellmed T-600-IUDT/AN, programme 2 for small parts; code B).

It is not allowed to use agents / disinfectants, which contain the following substances, as these may cause changes in the material:

- High-concentration organic and inorganic acids
- Chlorinated hydrocarbons
- 2-ethoxyethanol

When cleaning / disinfecting, the following procedures must be followed:

	Process	Time (sec.)
Zone 1	Pre-rinse, external, cold, 10 – 15°C Washing, acidic, external 35°C Draining time Re-rinse, external approx. 80°C Draining time Re-rinse, external approx. 80°C Draining time	45 120 10 *10 *15 *15 15
Zone 2	Washing, alkaline, external, 93°C Draining time Re-rinse, external, acidic, 90°C Draining time Re-rinse, external 90°C Draining time	135 10 10 15 15 15
Zone 3	Drying, external 100 – 120°C	200
Zone 4	Drying, external 100 – 120°C	200
	Door open / close & transport (sluice discharge)	60
	Cycle time overall ca.	290 ≈ 5 minutes

* When occupying the disinfection zone (washing zone 2), the re-rinse and draining times will depend on the respective objects being washed therein!

Sterilisation

Only previously cleaned and disinfected handles may be sterilised.

The handles are placed in a suitable sterilisation pack (one-way sterilisation pack, e.g. foil / paper sterilisation bags, single or double pack) in accordance with DIN EN 868 / ISO 11607 for steam sterilisation and then sterilised.

Use only the sterilisation procedure listed below for sterilisation. Other sterilisation procedures (e.g. ethylene oxide, formaldehyde and low-temperature plasma sterilisation) are not permissible.

Steam sterilisation procedure

Validated in accordance with DIN EN 554 / ISO 11134

Maximum sterilisation temperature 134°C

Proof in principle of the handles' suitability for effective sterilisation was provided using a fractional vacuum process (Euroselectomat 666 by MMM Münchner Medizin Mechanik GmbH, sterilising temperature 134°C, holding time 7 min.).

When applying other sterilization procedures, the suitability and effectiveness of the procedure must be proved in the context of the validation.

Inspection / durability

The handles should be inspected for damage and changed before re-use, if required.

The sterilisable handle sleeve must be disposed after 1000 sterilisation cycles or at the latest after 2 years and replaced with a new one.

The year of manufacture can be determined with the help of a stamping on the inner side of the handle sleeve (like shown in the photo). The stamping in the photo shows the number 12, which stands exemplary for the year 2012.

4.2 Lamp housing, protective lens and support system

The YLED-1F has a high-quality surface, which can be cleaned with conventional cleaning agents.

The clear protective cover 7 is made of a high-quality plastic. Pay attention to the following during cleaning:

- Never wipe over the clear protective cover 7 with a dry cloth (always clean with a wet cloth).
- Wipe the protective cover 7 after cleaning with an antistatic, non-fluffy cloth.

5. Maintenance

Preventive maintenance of the YLED-1F must be done every two years. This includes a technical and mechanical check-up.

Please observe also the mounting instructions and instructions for use of the carrying systems. These instructions can contain statements for different maintenance intervals.

Attention: Before dismantling the lamp, set the height adjustment of the spring arm to horizontal position (see mounting instructions „Portegra2“).

Attention: During all maintenance work the light must be disconnected from mains and secured against resetting.

5.1 Maintenance work

The following maintenance work / tests has / have to be done:

- check on defects in paint work;
- check on fissures at plastic parts;
- check on deformation of the suspension;
- check on loosened parts;
- check the connection between light and carrying systems;
- check and grease the securing segment;
- check the faultless function of the light;
- perform the electrical safety tests

For adjustments at the ceiling attachment please observe also the mounting instructions „Portegra2“.

6. Data

6.1 Technical data

	YLED-1F
Central light intensity at a distance of 1 meter	70,000 Lux
Light field diameter d_{10}	164 mm
Light field diameter d_{50}	87 mm
Light intensity with one shadower	0 %
Light intensity with two shadowers	47 %
Light intensity on the ground of a normed tube	100 %
Light intensity on the ground of a normed tube with one shadower	0 %
Light intensity on the ground of a normed tube with two shadowers	47 %
Illumination depth 20%	1750 mm
Illumination depth 60%	610 mm
Colour rendering index R_a	95
Colour rendering index R_9	90
Focusable light field size	14 - 25 cm
Colour temperature (Kelvin)	4100 K
Temperature increase in head area	0.5 °C
Electronic light intensity control at the light head	50 – 100 %
Number of LED's	17
Working distance	70 - 140 cm
Protection type of the lamp head per IEC 60529	IP 44
Dimension of light head	28 x 36 cm
Height adjustment	118 cm

Remark:

The technical data are subject to fluctuations. Due to manufacturing reasons the real values can slightly differ from the data mentioned above.

The values for R_a and R_9 can differ with approx $\pm 5\%$.

The values for the colour temperature can differ with approx $\pm 200K$.

6.2 Electrical data

YLED-1F	
Power consumption	0,5 A
Operating voltage AC	100 - 240 V

6.3 Information regarding the electrical installation

When turned ON, the examination light YLED-1F is exposed to a current peak.

The electronic power supply is integrated in the light head.

The examination light YLED-1F can be connected directly to the mains. The wide range input can be operated with an input voltage of 100 – 240 VAC, 50 – 60 Hz.

Warning!

The light is class I equipment. In order to avoid the risk of an electric shock, the equipment must be connected to a mains supply with protective earth.

6.4 Weights

Light	Weight
YLED-1F	2,4 kg

6.5 Environmental conditions

Operation



	Min.	Max.
Temperature	+10°C	+30°C*
Relative atmospheric humidity	30 %	75 %
Air pressure	700 hPa	1060 hPa

* in case of higher temperatures please contact us

Transport / storage

	Min.	Max.
Temperature	-10°C	+50°C
Relative atmospheric humidity	20 %	90 %
Air pressure	700 hPa	1060 hPa

References on the package

6.6 General remarks

When using more than one examination light at the same time (light combinations), please make sure that the total radiation intensity in field is less than 1000 W/m^2 , to avoid further temperature increase in the wound field.

When installing an examination light, its fail-safety must be guaranteed according to DIN VDE 0100-710 (former DIN VDE 0107).

The test certificate for the electrical safety test can be requested when needed. Please provide the serial number of the respective light.

In case of a collective wiring of further lights or devices at installation, chapter 16 of the European standard EN 60601-1:2013 must be applied and eventually it has to be checked if the requirements are met.

The light must be tested according to EN 62353 at commissioning.

7. CE-mark

The products YLED-1F comply with the European Directive 93/42/EEC for medical devices of the European Community's Council.
Dr. Mach applies the standard EN 60601-2-41.
Dr. Mach GmbH is certified according to DIN EN ISO 13485:2016 + AC: 2016.

8. Manufacturer

Dr. Mach GmbH & Co. KG, 85560 Ebersberg - Germany

9. Disposal

The examination light does not contain any dangerous goods.
The components of the examination light should be properly disposed at the end of its shelf-life.
Make sure, that the materials are carefully separated.
The electrical conducting boards should be submitted to an appropriate recycling proceeding.
The rest of the components should be disposed according to the contained materials.

10. Electromagnetic compatibility

The Dr. Mach OT- and examination lights are subject to special preventive measures regarding the Electromagnetic compatibility and must be installed according to the EMC-instructions mentioned in the accompanying documents.

The function of the OT- and examination lights can be affected by portable and mobile HF-communication or other HF devices.

If the essential performance should be lost or degraded due to electromagnetic disturbances it is to be expected that the provisioning of the illumination, for a short time, is not ensured.

To ensure the basic safety and the essential performance of the YLED-1F regarding electromagnetic disturbances during the expected service life, the maintenance has to be performed within the defined interval and according to the instructions in the technical service manual. Only spare parts specified by the legal manufacturer may be installed.

The Examination light YLED-1F has been tested under professional healthcare facility environment.

The use of other equipment than the equipment mentioned in chapter 2. leads to an increased emission or to a reduced interference resistance of the device.

For the intended use of the examination light YLED-1F it is required that the YLED-1F is not mounted immediate and near to other devices or mounted together with other devices. If operating the light is obligatory near other devices or together with other devices, the functions of the light YLED-1F must be observed.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the YLED-1F, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Table 2 – Emission limits per environment

Phenomenon	Professional healthcare facility environment ^{a)}	Home healthcare environment ^{a)}
Conducted and radiated RF emissions	Class B; Group 1 (according CISPR 11)	CISPR 11 ^{c), d)}
Harmonic distortion	See IEC 61000-3-2 ^{b)}	See IEC 61000-3-2
Voltage fluctuations and flicker	See IEC 61000-3-3 ^{b)}	See IEC 61000-3-3
^{a)} See 8.9 of IEC 60601-1-2:2014 for information about the environments of intended use. ^{b)} This test is not applicable in this environment unless the ME EQUIPMENT and ME SYSTEMS used there will be connected to the public mains network and the power input is otherwise within the scope of the Basic EMC standard. ^{c)} ME EQUIPMENT and ME SYSTEMS intended for use in aircraft shall meet the RF emissions requirements of ISO 7137. The conducted RF emissions test is applicable only to ME EQUIPMENT and ME systems that are intended to be connected to aircraft power. ISO 7137 is identical to RTCA DO-160C:1989 and EUROCAE ED-14C:1989. The latest editions are RTCA DO-160G:2010 and EUROCAE ED-14G:2011. Therefore, use of Section 21 (and category M) of a more recent edition, e.g. [39] or [40], should be considered. ^{d)} Standards applicable to other modes or EM environments of transportation for which use is intended shall apply. Examples of standards that might be applicable include CISPR 25 and ISO 7637-2.		

Table 4 – Enclosure port

Phenomenon	Basic EMC standard or test method	Immunity test levels	
		Professional healthcare facility environment	Home healthcare environment
Electrostatic Discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV , ± 15 kV air	
Radiated RF EM fields ^{a)}	IEC 61000-4-3	3 V/m ^{f)} 80 MHz - 2,7 GHz ^{b)} 80 % AM at 1 kHz ^{c)}	10 V/m ^{f)} 80 MHz - 2,7 GHz ^{b)} 80 % AM at 1 kHz ^{c)}
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See Table 9	
Rated power frequency magnetic fields ^{d)e)}	IEC 61000-4-8	30 A/m ^{g)} 50 Hz or 60 Hz	

^{a)} The interface between the patient physiological signal simulation, if used, and the ME EQUIPMENT and ME SYSTEM shall be located within 0,1m of the vertical plane of the uniform field area in one orientation of the ME EQUIPMENT and ME SYSTEMS.

^{b)} ME EQUIPMENT and ME SYSTEMS that intentionally receive RF electromagnetic energy for the purpose of their operation shall be tested at the frequency of reception. Testing may be performed at other modulation frequencies identified by the risk management process. This test assesses the basic safety and essential performance of an intentional receiver when an ambient signal is in the passband. It is understood that the receiver might not achieve normal reception during the test.

^{c)} Testing may be performed at other modulation frequencies identified by the risk management process.

^{d)} Applies only to ME EQUIPMENT and ME SYSTEMS with magnetically sensitive components or circuitry.

^{e)} During the test, the ME EQUIPMENT or ME SYSTEMS may be powered at any nominal input voltage, but with the same frequency as the test signal (see Table 1).

^{f)} Before modulation is applied.

^{g)} This test level assumes a minimum distance between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic field of at least 15 cm. If the risk analysis shows that the ME EQUIPMENT or ME SYSTEM will be used closer than 15 cm to sources of power frequency magnetic field, the immunity test level shall be adjusted as appropriate for the minimum expected distance.

Table 5 – Input a.c power PORT (1 of 2)

Phenomenon	Basic EMC standard	Immunity test levels	
		Professional healthcare facility environment	Home healthcare environment
Electrostatic fast transients bursts ^{a) l) o)}	IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	
Surges ^{a) b) j) k) o)} Line – to- line	IEC 61000-4-5	± 0,5 kV, ± 1 kV	
Surges ^{a) b) j) o)} Line – to- ground	IEC 61000-4-5	± 0,5 kV, ± 1 kV, ± 2 kV	
Conducted disturbances induced by RF fields ^{c) d) o)}	IEC 61000-4-6	3 V ^{m)} 0,15 MHz - 80 MHz 6 V ^{m)} in ISM bands between 0,15 MHz and 80MHz ⁿ⁾ 80 % AM at 1 kHz ^{e)}	3 V ^{m)} 0,15 MHz - 80 MHz 6 V ^{m)} in ISM and amateur radio bands between 0,15 MHz and 80MHz ⁿ⁾ 80 % AM at 1 kHz ^{c)}
Voltage dips ^{f) p) r)}	IEC 61000-4-11	0 % U_T ; 0,5 cycle ^{g)} At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° ^{q)} 0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycles ^{h)} Single phase: at 0°	
Voltage interruptions ^{f) i) o) r)}	IEC 61000-4-11	0 % U_T ; 25/30 cycles ^{h)}	

- a) The test may be performed at any one power input voltage within the ME EQUIPMENT OR ME SYSTEM RATED voltage range. If the ME EQUIPMENT OR ME SYSTEM is tested at one power input voltage, it is not necessary to re-test at additional voltages.
- b) All ME EQUIPMENT and ME SYSTEM cables are attached during the test.
- c) Calibration for current injection clamps shall be performed in a 150.Ω system.
- d) If the frequency stepping skips over an ISM or amateur band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- e) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- f) ME EQUIPMENT and ME SYSTEMS with a d.c. power input intended for use with a.c.-to-d.c. converters shall be tested using a converter that meets the specifications of the MANUFACTURER of the ME EQUIPMENT or ME SYSTEM. The IMMUNITY TEST LEVELS are applied to the a.c. power input of the converter.
- g) Applicable only to ME EQUIPMENT and ME SYSTEMS connected to single-phase a.c. mains.
- h) E.g. 10/12 means 10 periods at 50 Hz or 12 periods at 60 Hz.
- i) ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase shall be interrupted once for 250/300 cycles at any angle and at all phases at the same time (if applicable). ME EQUIPMENT and ME SYSTEMS with battery backup shall resume line power operation after the test. For ME EQUIPMENT and ME SYSTEMS with RATED input current not exceeding 16 A, all phases shall be interrupted simultaneously.
- j) ME EQUIPMENT and ME SYSTEMS that do not have a surge protection device in the primary power circuit may be tested only at ± 2 kV line(s) to earth and ± 1 kV line(s) to line(s).
- k) Not applicable to CLASS II ME EQUIPMENT and ME SYSTEMS.

Table 5 (2 of 2)

- l) Direct coupling shall be used..
- m) r.m.s., before modulation is applied.
- n) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50 MHz to 54,0 MHz.
- o) Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase and ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase.
- p) Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase.
- q) At some phase angles, applying this test to ME EQUIPMENT with transformer mains power input might cause an overcurrent protection device to open. This can occur due the magnetic flux saturation of the transformer core after the voltage dip. If this occurs, the ME EQUIPMENT or ME SYSTEM shall provide BASIC SAFTEY during and after the test.
- r) For ME EQUIPMENT and ME SYSTEMS that have multiple voltage settings or auto ranging voltage capability, the test shall be performed at the minimum and maximum RATED input voltage. ME EQUIPEMT and ME SYSTEMS with a RATED input voltage range of less than 25 % of the highest RATED input voltage shall be tested at one RATED input within the range. See table 1 Note c) for examples calculations.

Table 9 – Test specifications for enclosure port immunity to RF wireless communications equipment

Test Frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum Power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 - 390	TETRA 400	Pulse modulation ^{b)}	1,8	0,3	27
450	430 - 470	GMRS 460, FRS 460	18 Hz FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0,3	28
710	704 – 787	LTE Band 13, 17	Pulse modulation ^{b)}	0,2	0,3	9
745						
780			217 Hz			
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)}	2	0,3	28
870						
930			18 Hz			
1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)}	2	0,3	28
1845						
1970			217 Hz			
2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28
5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation ^{b)}	0,2	0,3	9
5500						
5785			217 Hz			

NOTE If necessary to achieve the immunity test level, the distance between the transmitting antenna and the ME equipment or ME system may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

^{a)} For some services, only the uplink frequencies are included.

^{b)} The carrier shall be modulated using a 50 % duty cycle square wave signal.

^{c)} As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

11. Overview of the light functions

Switching the examination light ON/OFF is done by...

...pressing
the ON/OFF button

Focusing

...by turning
the handle

The light intensity control is adjusted by...

...pressing
the arrow buttons