YLED-1F

DE	deutsch Gebrauchsanweisung	FI	suomi Käyttöohjeet	SE	svenska Bruksanvisning
EN	english User manual	HU	magyar Használati utasítás	SK	slovenčina Návod na použitie
FR	français Mode d'emploi	HR	hrvatski Uputa za uporabu	SL	slovenščina Navodila za uporabo
IT	italiano Istruzioni per l'uso	LT	lietuvių Naudojimo instrukcijos	TR	Türkçe Kullanım Kılavuzu
ES	español Manual de instrucciones	LV	latviešu Lietošanas instrukcija	RU	Русский Инструкция по эксплуатации
BG	български език Инструкция за употреба	NL	nederlands Gebruikershandleiding	ID	bahasa Indonesia Petunjuk penggunaan
CS	Česky Návod k použití	NO	Norsk Bruksanvisning	JP	日本語 取扱説明書
DA	dansk Brugsanvisning	PL	wersja polska Instrukcja obsługi	KO	한국어 사용 설명서
EL	Ελληνικά Οδηγίες χρήσης	PT	português Manual de instruções	VI	Tiếng Việt Hướng dẫn về cách sử dụng
ET	eesti Kasutusjuhend	RO	română Manual de utilizare	ZH	中文 使用说明书

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Congratulations on acquiring a new YLED-1F examination light.

Please read these Instructions for use very carefully.

1. Instructions for safe use

1.1 Intended user

The YLED-1F is a Class I medical product and may only be operated by trained medical staff.

1.2 Information and obligation of the user to check the product

Pay attention to the instructions for use when handling the lamp. These Instructions for use are part of the product and must therefore be stored in a place close to the product in order for the safety instructions and important information to be consulted at any time.

Make sure that the lamp is in satisfactory working order before every use. If there is obvious damage, unusual operating conditions, etc., the lamp must not be used.

1.3 Availability of the instructions for use

These instructions for use and a detailed manual with further tips on how to use the light can be found online at the following link:

<https://mavig.com/eifu/>

1.4 Intended use / contra-indications

The YLED-1F examination light is designed to illuminate an examination site in medical facilities (e.g. in a laboratory, in hospitals or doctor's practice) with focused, low-glare, shadow-free light. It enables the user to perform a diagnosis or carry out medical interventions. The YLED-1F light is an operating light that is not fail-safe when used as a single light. This device is not designed to operate in areas which are subject to explosions.

Continuous illumination of the open human eye should be avoided when illuminating the face area.

1.5 Technical data

Class of protection	I
IP protection class	IP 44
Input voltage (light body) ¹	100-240 V AC, 50/60 Hz
Current	0.75-0.5 A max.
Examination time	Continuous operation possible
Expected life ²	10 years

¹ The electronic power supply unit is integrated in the light head.

² At the end of the expected (designed) service life, the lamp must be serviced more frequently for safe operation (see the manual for details).

1.6 Lighting technical data

	YLED-1FF
Central illuminance (distance 1 m)	70,000 lux
Light field diameter d10	164 mm
Light field diameter d50	87 mm
Residual light intensity (one shadower)	0 %
Residual light intensity (two shadowers)	47 %
Residual light intensity (standardized tube)	100 %
Residual light density (standardized tube, one shadower)	0 %
Residual light density (standardized tube, two shadowers)	47 %
Depth of illumination (20%)	1750 mm
Depth of illumination (60%)	610 mm
Total irradiance in the field (distance 1 m)	250 W/m ²
Total irradiance in the field (distance 0.90 m)	268 W/m ²

For a complete overview of technical and lighting specifications, refer to the manual.

1.7 Installation/Maintenance/Repair

The lamp may only be installed, maintained or repaired by the manufacturer or by specially trained staff. Maintenance must be carried out at least every two years

1.8 Environmental conditions for operation

Ambient temperature:	+5°C to +40°C
Relative air humidity	30% to 75% RH
Air pressure	700 hPa to 1060 hPa

1.9 Reporting obligations

All serious incidents which have occurred in connection with the product must be reported to the manufacturer and the competent authority.

2. Images on the device

This symbol indicates that the instructions for use must be followed.

This symbol indicates that documents are available digitally

Serial number of the product

Part number of the product

Address of the manufacturer

Date of manufacture and country of manufacture

EC conformity symbol

This symbol indicates that this is a medical device

Unique device identifier (ID) of the product

NRTL Test mark
The lamp was tested by a 'Nationally Recognized Testing Laboratory'

Reference to China RoHS / pollution control logo

Instruction for disposal of the device

3. Safety instructions

	The instructions for use must be observed for safe handling of the lamp
	To avoid the risk of an electric shock, this device must only be connected to a supply network that has a protective earthing conductor.
	A primary-side ON/OFF switch to isolate the system from the supply network must be provided on site. The switch must meet the requirements of IEC 61058-1 for nominal voltage peaks of 4 kV.
	This device is not designed for operation in environments enriched with oxygen.
	The lamp may only be used for the intended purpose. Otherwise, the manufacturer will not be liable for personal injury or damage to property.
	The lamp is equipped with a sterilisable handle at the factory and must only be used with this handle.
	Changes to the light are prohibited and will invalidate the manufacturer's certificate of conformity and all warranty claims.
	Installation, maintenance and repair work may only be carried out by the manufacturer or by specially trained staff.
	Maintenance must be carried out on the light at least every two years

Additional equipment that is connected to medical electrical equipment, must conform to the relevant IEC- or ISO standards (e.g. IEC 60950 or IEC 62368 for data processing equipment). Moreover, all configurations must meet the requirements for medical electrical systems (see Section 16 of the latest version of IEC 60601-1). Anyone who connects additional equipment to medical electrical equipment, is configuring a medical system and is therefore responsible for ensuring that the system meets the requirements for medical electrical systems. In case of doubt, contact your local representative or our technical customer services.

The simultaneous use of several lights to illuminate a wound area may result in the maximum allowed energy input being exceeded (1,000 W/m²) and thus excessive heat development. It is the user's responsibility to ensure that the maximum allowed limit is not exceeded.

The unprotected human eye can be damaged by direct light. Do not look directly into the light beam of the lamp. Do not point the light beam at the patient's unprotected eye continuously.

When positioning the light body, there is a risk of injury (e.g. crushing) and collisions with other objects (inventory) or walls.

Parts that fall off could injure the patient or lead to an infection of the wound area.

Do not remove the rating plate or the warning labels.

It is forbidden to carry out servicing or repair activities whilst the lamp is in use.

4. Operating the YLED-1F light


ON/OFF switch of the light (to switch off, hold for one second) (1)

Adjust the brightness in six levels between 50 and 100% (2 and 3) (lowest level = endoscopy light)

Display of the set brightness (4), the uppermost LED will flash in Standby

If only the top LED lights up, the light is in endoscopy mode

Focus the light field by turning the handle (5)

Position the lamp body using the handle (5)

5. Cleaning and disinfecting

Cleaning and disinfecting work must only be done by trained staff. The respective requirements must be observed for all cleaning and disinfection work (details can be found in the manual).

Housing/protective screen

The housing and the protective screen of the lamp body can be cleaned and disinfected with many common/commercially available materials. **Do not use cleaning agents or disinfectants containing active substances based on biguanides, phenols, chlorides or percompounds.**

Furthermore, only cleaning agents approved for polycarbonate (PC) may be used to clean the protective screen. To protect against mechanical damage, always use a damp cloth (never a dry one) to wipe the protective screen and after cleaning, wipe with an anti-static agent (lint-free cloth).

Sterilisable handle

The handle must be cleaned/disinfected before each use. It can be steam sterilised (max. 200 sterilisation cycles for max. 5 minutes at a max. temperature of 134°C); details can be found in the manual.

Before installing the handle, check it for visible damage, fouling and the specified manufacturing date. Do not use damaged or dirty handles or handles that are more than two years old.

6. Faults

In the event of unusual examination conditions, the lamp must not be used, as safe operation cannot be guaranteed. For troubleshooting, disconnect the lamp from the mains for about 30 seconds. In the event of continuous faults, a suitably trained service technician must be contacted.

7. Information on electromagnetic compatibility

Medical electrical equipment is subject to special precautionary measures with regard to electromagnetic compatibility (EMC). It may only be installed and put into service in accordance with the EMC instructions in the accompanying documents. The YLED-1F operating light has been tested for use in professional equipment of the health care system.

	Portable and mobile RF communications equipment can affect medical electrical equipment and must not be used within 30 cm of the light, including the cable.
	The use of this equipment immediately next to other equipment or with other equipment in stacked form should be avoided since this may result in faulty operation. Should use in the aforementioned manner nevertheless be necessary, this device and the other equipment should be kept under observation to ensure that they are working properly.
	The use of different accessories, converters or cables to those that the manufacturer of this device has stipulated or made available may result in increased electromagnetic interference or reduced immunity to electromagnetic interference and to faulty operation.
	Furthermore, the light must not be operated if the housing, cable or electromagnetic shielding is damaged.

Additional information on electromagnetic compatibility can be found in the manual.

8. Disposal

The light does not contain any harmful substances. The components of the light should be disposed of appropriately at the end of the product's life.

Take care that the material is carefully separated: The electrical circuit boards should be recycled appropriately. The housing of the light and the other components should be disposed of according to the materials they contain.

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