Videosysteme

DE	deutsch Gebrauchsanweisung	FI	suomi Käyttöohjeet	SE	svenska Bruksanvisning
EN	english User manual	HU	magyar Használati utasítás	SK	slovenčina Návod na použitie
FR	français Mode d'emploi	HR	hrvatski Uputa za uporabu	SL	slovenščina Navodila za uporabo
IT	italiano Istruzioni per l'uso	LT	lietuvių Naudojimo instrukcijos		
ES	español Manual de instrucciones	LV	latviešu Lietošanas instrukcija		
BG	български език Инструкция за употреба	NL	nederlands Gebruikershandleiding		
CS	Česky Návod k použití	NO	Norsk Bruksanvisning		
DA	dansk Brugsanvisning	PL	wersja polska Instrukcja obsługi		
EL	Ελληνικά Οδηγίες χρήσης	PT	português Manual de instruções		
ET	eesti Kasutusjuhend	RO	română Manual de utilizare		

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Congratulations on your purchase of a new video system

Please read these instructions for use very carefully.

1. Instructions for safe use

1.1 Intended user

The video system is a Class I medical product and may only be operated by trained medical staff.

1.2 Information obligation and user obligation to check the product

To operate the video system, these instructions for use must be observed. These instructions for use form part of the product and must therefore be stored in a place close to the product so that the safety instructions and important information can be consulted at any time.

Before each use, make sure that the video system is in perfect condition. In the event of obvious damage, unusual operating conditions, etc., the system must not be used!

1.3 Availability of the instructions for use

These instructions for use and a detailed manual with further tips on how to use the video system optimally can be found online at the following link:

<https://dr-mach.de/login/machvideo.html>

1.4 Intended use / contra-indications

This video system is designed to record a video image from an examination or operating theatre site (e.g. in a laboratory, a hospital or a doctor's practice) and to make it available as a video signal via a standardised interface. The video signal can be used for medical education and training or for documentation.

It is not permitted to use the video image to support a treatment or to make a diagnosis.

1.5 Technical data

Video system	HDMV	HDMV F	HDMV F 4K
Video resolution	HD 1080i/60, 1080i/50	HD 1080p/60, 1080p/50	4K 2160p/30, 2160p/25
Zoom range	30x optical, 12x digital	30x optical, 12x digital	30x optical, 12x digital
Video output	1 x HD - SDI	1 x HDMI	1 x HDMI
Power supply / power consumption	24-30 V DC / 10.5W	24-30 V DC / 11.5W	24-30 V DC / 13W
Expected life ⁽¹⁾	10 years	10 years	10 years

¹After the intended (designed) service life has expired, the video system must be serviced more frequently to ensure safe operation (see manual for details).

1.6 Installation/Maintenance/Repair

The video system may only be installed, maintained or repaired by the manufacturer or by specially trained staff. Maintenance must be carried out at least every two years

1.7 Environmental conditions for operation

Ambient temperature: +5 °C to +40 °C
Relative air humidity: 30% to 75% RH
Air pressure: 700 hPa to 1060 hPa

1.8 Reporting obligations

All serious incidents that occur in connection with this product must be reported to the manufacturer and the competent authority.

2. Images on the device

This symbol indicates that the instructions for use must be followed.

Product serial number

Product item number

Manufacturer's address

Date of manufacture and country of manufacture

EC conformity symbol

This symbol indicates that this is a medical device

Unique device identifier (UDI) of the product

Reference to China RoHS / pollution control logo

Information about device disposal


3. Safety instructions

	This symbol indicates possible sources of danger. Please also note the safety instructions and the hazard specification in the associated installation and operating instructions for the support arm system.
	To safely operate the video system, these instructions for use must be observed.
	To avoid the risk of an electric shock, this device must only be connected to a supply network that has a protective earthing conductor.
	A primary-side ON/OFF switch to isolate the system from the supply network must be provided on site. The switch must meet the requirements of IEC 61058-1 for nominal voltage peaks of 4kV.

	This device is not designed for operation in environments enriched with oxygen.
	The video system may only be used in accordance with its intended purpose. Otherwise, the manufacturer will not be liable for personal injury or damage to property.
	The external camera holder is factory-equipped with a sterilisable handle and may only be used with this.
	Use only the mains units approved or supplied by the manufacturer. Failure to do so will invalidate the conformity of the product and all warranty claims.
	Modifications to the video system are prohibited and will void the manufacturer's certificate of conformity and all guarantee/warranty claims.
	Installation, maintenance and repair work may only be carried out by the manufacturer or by specially trained staff.
	Maintenance of the video system must be carried out at least every two years.
	Additional equipment that is connected to medical electrical equipment must conform to the relevant IEC- or ISO standards (e.g. IEC 60950 or IEC 62368 for data processing equipment). Moreover, all configurations must meet the requirements for medical electrical systems (see Section 16 of the latest version of IEC 60601-1). Anyone who connects additional equipment to medical electrical equipment is configuring a medical system and is therefore responsible for ensuring that the system meets the requirements for medical electrical systems. In case of doubt, contact your local representative or our technical customer service department.
	When positioning the camera, there is a risk of injury (e.g. by crushing) and collisions with other objects (inventory) or walls.


	Parts that fall off could injure the patient or lead to an infection of the wound area.
	Do not remove the rating plate or the warning labels.
	It is forbidden to carry out maintenance or repair activities whilst the video system is in use.
	Simultaneous touching of parts on the luminaire and the patient is not permitted.

4. Operating the video system

Installing the camera in a light:

Remove cover by pressing the locking pins on the side.


Installation:

Push the camera into the opening until both locking pins engage securely

Check that the camera is firmly seated!

Removal:

Press the locking pins on the side and carefully remove the camera

Installing the camera in the external camera holder:

Slide camera into camera holder
Observe alignment of plug 1 and socket!

Secure the camera by tightening the pin 2

Remove camera by loosening pin 2

Check that the camera is firmly seated!

Camera controls:

ON/OFF:

Switch between operation and stand-by

Rotate:

Rotate video image

Zoom:

Magnification of the video image

White balance:



adjusting the colour temperature to the actual environment

Still image

Focusing:


Autofocus – Switch on/off

Near/far – Adjust to near or far objects


Exposure:
 AE – Switching the automatic exposure control on/off.
 Sun/Moon – Increase or decrease the brightness

Connection of the radio receiver:
 Pull off the top cover **3**

Insert the receiver module **4** into the wall bracket

Pairing (HDMV F and HDMV F 4K only):
 Insert the camera into the side connector while supporting it

Start pairing by pressing the red pairing button

After approx. three minutes, switch on the camera via the remote control (the image should then be displayed on the monitor).

5. Cleaning and disinfecting

Cleaning and disinfecting work must only be done by trained staff. The respective requirements must be observed for all cleaning and disinfection work (details can be found in the manual).

Camera

The video system can be cleaned and disinfected using many common / commercially available agents. **Cleaning agents or disinfectants with active ingredients based on chloromethylphenols, phenols or sodium hydroxide must not be used!**

Only wipe the camera unit with a damp cloth (not dripping).

Furthermore, only cleaning agents approved for polycarbonate (PC) may be used to clean the lens. To protect against mechanical damage, always use a damp cloth (never a dry one) to wipe the lens and, after cleaning, wipe with an anti-static agent (lint-free cloth).

Sterilisable handle

The handle must be cleaned/disinfected before each use. It can be steam sterilised (max. 200 sterilisation cycles for max. five minutes at a max. temperature of 134°C; details can be found in the manual).


Before installing the handle, check it for visible damage, soiling and the specified manufacturing date. Do not use damaged or dirty handles or handles that are more than two years old.

6. Faults

In the event of unusual operating conditions, the video system must not be used, as safe operation cannot be guaranteed. For troubleshooting, disconnect the video system from the mains for approx. 30 seconds; if errors persist, contact a suitably trained service technician.

7. Notes on electromagnetic compatibility

Medical electrical equipment is subject to special precautionary measures with regard to electromagnetic compatibility (EMC). It may only be installed and put into service in accordance with the EMC instructions in the accompanying documents. The video system has been tested for use in professional healthcare facilities.

	The use of this equipment immediately next to other equipment or with other equipment in stacked form should be avoided since this may result in faulty operation. Should use in the aforementioned manner nevertheless be necessary, this device and the other equipment should be kept under observation to ensure that they are working properly.
	The use of different accessories, converters or cables to those that the manufacturer of this device has stipulated or made available may result in increased electromagnetic interference or reduced immunity to electromagnetic interference and to faulty operation.
	Furthermore, the video system must not be operated if the housing, cable or electromagnetic shielding is damaged.

Additional information on electromagnetic compatibility can be found in the manual.

8. Disposal

The video system does not contain any harmful substances. At the end of the product's life, its components should be disposed of properly.

Take care that the materials are carefully separated: The electrical circuit boards should be recycled appropriately. The video system's remaining components should be disposed of according to the materials.

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<https://dr-mach.de/login/machvideo.html>

1.4 Usage conforme/contre-indications

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