


EN

**OT50001**  
**OT90001**

**MAVIG**



**STEA-OT5** (accessory) Disposable, sterile cover for the whole shield

**Article Codes:** OT50001  
OT90001

Shield complete with Portegra2 Suspension Arm (extension arm length 75 cm)  
Shield complete with Portegra2 Suspension Arm (extension arm length 95 cm)

IEC 60601-1

Made in Germany




## OT50001/OT90001 – X-Ray Protective Shield

MAVIG's products are designed to offer solid protection and optimized functionality. The high quality acrylic shield allows for an excellent transmittance and unimpeded view. All edges and corners have been rounded and shatter-proof materials are used to minimize risk of injury.

MAVIG's patented Portegra2 suspension arm provides its industry known safety features and movement flexibility. The centrally guided shield with a special patient cut-out can be easily positioned over the patient and quickly removed in case of an emergency.

- **Lead Acrylic Shield** W x H: 60 x 76 cm (23.6 x 29.9 in)  
Weight: 9.7 kg (21.4 lbs)
- **Lead Equivalent** 0.50 mm Pb
- **Suspension Arm** Extension / Spring Arm  
Lengths: OT50001 75 cm / 91 cm  
OT90001 95 cm / 91 cm  
Colour: RAL 9002
- **Replacement Parts** E-OT25B05: X-ray protective shield with connection bar (500 mm)



### MAVIG GmbH

X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone +49 (0) 89 / 420 96 - 0  
Fax +49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com