

EN

BRYGGA®

MAVIG



Model: BRYGGA Carbon brackets, weight relief for the shoulder when wearing X-ray protective aprons

Made in Sweden

BRYGGA® – Shoulder Weight Relief for Apron Users

Wearing the BRYGGA shoulder bridge underneath the X-ray apron is an easy-to-use and very efficient method of reducing physical strain and muscular fatigue. This patented solution from Sweden is a unique way of relieving pressure normally caused by the weight of the apron loaded onto the shoulder area.

Spreading the weight more evenly and resting on the bone structure rather than the muscles and nerve-centre, the use of BRYGGA is particularly beneficial during long procedures.

- **Patented solution**
- **Recommended for use with MAVIG apron models RA631 and RA632**
- **Does not hinder the user or interfere with procedures**
- **Promotes an ergonomically correct posture**
- **Greatly reduces the risk of tension headache and numbness of the arms**
- **Full freedom of movement of the neck and arms**
- **Easy cleaning** (just wipe off with any cleaning agent)
- **Available in sizes S/M and L**
- **Simple adjustment**
- **Very comfortable**
- **Outer material 100% durable carbon fiber**
- **Padding foam:** **Outside** - Evazote, density 50
Inside - Polyethylene, density 70

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com