

EN

BR331

Design Blue-transparent

MAVIG



Model: BR331

Glasses with front X-ray protection of 0.75 mmPb and side X-ray protection of 0.50 mmPb, universal size, weight approx. 85 g

Made in Germany

CE 0493

BR331 – X-Ray Protective Glasses

The X-ray protective glasses have flat glasses and a radiation protection of 0.75 mmPb. The special plastic shells, which are lined on the inside with coated X-ray protective material up to the level of the earpieces, offer lateral protection of 0.50 mmPb. The frame is made of an extremely durable plastic which ensures easy cleaning and good resistance to chemicals (disinfectant agents). The BR331 comes with a storage case and a soft head band.

MAVIG uses only protective lenses of the highest quality. This guarantees the best light transmission values and a maximum lens break resistance. An ideal composition of the protective material results in a high level X-ray absorption, and optimized production processes achieve an absolute homogeneous allocation of the glass additives.

The quality features of our glasses make them perfect for the use in the medical environment:

- Great frame stability through high-grade synthetic materials
- Optimized balance of sturdy, light-weight frame and lenses
- Minimized risk of injury due to shape and used components
- Low allergy risk due to carefully chosen materials
- Easy cleaning and superior resistance against chemical influences (disinfection)
- Wide and flexible head band for a comfortable fit
- For this model prescription lenses are available: single, bifocal, and progressive
- Conformity: DIN EN 61331-1:2016 / DIN EN 61331-3:2016

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com