



Model: RP2705 Breast protection for female patients during mammographic radiation procedures

Made in Germany


CE 0302

## RP270 – Breast Protection for Women

The mammary gland shield efficiently protects your patients during mammograms.

The shield consists of two equal parts, each 18 x 40 cm (W x H), connected by a lengthadjustable hook-and-loop fastener for comfortable donning and easy adjustment.

- **Protective Material** Novalite radiation protective material
- **Lead Equivalent** 0.50 mm Pb
- **Size** Universal
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Plus (design Happy only), available in the colours:  
*Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or designs: Happy*
- **Edging** All edges are covered with a seem band.



**MAVIG GmbH**  
X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone +49 (0) 89 / 420 96 - 0  
Fax +49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com