



Models:	RP2761E	Testicle Protection for men
	RP2761K	Testicle Protection for boys
	RP2761G6232	Testicle Shield Set for Children aged 1 to 10

Made in Germany

CE 0302

RP276 – Testicle Shield for Men/Boys/Children

This male gonad protection model covers the penis and scrotum all around with a lead equivalent of 1.00 mm Pb. Of course, MAVIG has made sure that there is no risk of injury to the patient when using the protection.

The hook-and-loop fastener provides a good hold when adjusted correctly. For cleaning and disinfection, this model series can be unfolded using the pressstuds. We also have suitable disposable hygiene bags* in our range.

- **Protective Material** Standard lead
- **Lead Equivalent** 1.00 mm Pb
- **Outer Cover/Colour** Soft PVC cover in light grey (boys/men);
Textile cover made of special, easy to clean material
in design Happy (infants/baby/children)
- **Sizes**
 - RP2761E - for men (dimensions approx. 155 x 120 x 90 mm)
 - RP2761K - for boys (dimensions approx. 110 x 90 x 65 mm)
 - RP2761G6232 - Set with 3 protections from infant to child sizes

Hint: Different hygienic cover sets are available for the sizes E and K, which each contain 100 covers. The set RP-HYG-CEK contains 100 covers for size E and 100 covers for size K.

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com