


EN

RP278

MAVIG



Models:	RP2781L	Ovarian Protection, Large	RP2781G	Ovarian Protection Set, all three sizes
	RP2781M	Ovarian Protection, Medium		
	RP2781S	Ovarian Protection, Small		

Made in Germany

CE 0302

RP278 – Ovarian X-Ray Protection for Female Supine Patients

The concept of this MAVIG ovarian shield is unsurpassed and meets all the requirements of radiological practice. Radiation exposure of the female gonads is significantly reduced without compromising the diagnostic technique.

For optimal fine adjustment, each model can be varied again by means of the modular principle. The RP278 is, of course, easy to clean and disinfect. A complete set* offers you the best possible protection for all your supine patients.

- **Protective Material** Standard lead
- **Lead Equivalent** 1.00 mm Pb
- **Outer Cover/Colour** Soft PVC cover in light grey
- **Sizes**
 - Large – Height 12 cm
 - Medium – Height 9 cm
 - Small – Height 7 cm

MAVIG GmbH

X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com

* Also available as the complete set (all 3 sizes) with price saving = RP2781G