

EN RP280

MAVIG



Models:	RP2801E	Testicle Protection for Adults (M)	RP-HYG-CM	Hygienic Covers (100 pcs, adults)
	RP2801K	Testicle Protection for Boys (K)	RP-HYG-CK	Hygienic Covers (100 pcs, boys)

Made in Germany

CE 0302

RP280 – Testicle X-Ray Protection for Men/Boys

The penis and scrotum are protected with a lead equivalent of 1.00 mm Pb. The padded, elastic pouch closure provides gentle pressure for a secure and correct fit throughout the examination. The rounded design poses no risk of injury to the patient.

It is easy to put on and the skin-friendly coating of the MAVIG RP280 testicle capsules offers seamless hygiene inside and out. Cleaning and disinfection are easy.

- **Protective Material** Lead rubber singleply material
- **Lead Equivalent** 1.00 mm Pb
- **Outer Cover/Colour** Skin friendly coating on a natural rubber base/
Medico-white
- **Sizes** RP2801E: Adults (size M, approx. 135 x 105 x 60 mm)
RP2801K: Boys (size K, approx. 105 x 85 x 55 mm)
- **Hygienic Covers** Different hygienic cover sets are available for the sizes M and K

Hint: Different hygienic cover sets are available for the sizes E and K, which each contain 100 covers. The set RP-HYG-CEK contains 100 covers for size E and 100 covers for size K.

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com