

Models:

RP6435SGonad apron (toddlers)**RP6435M**Gonad apron (children)

RP6435L RP6435X Gonad apron (adults) Gonad apron (X-large)

Made in Germany

C€ 0302

RP643 – Gonad Apron for Patients

An article for daily use makes special demands on material quality and workmanship but the functionality for user-friendly use must also be allowed for. Even with the high patient flow in today's radiology practices, the RP643 is ideally suited, as the X-ray protection is quick and practical to put on and take off. It couldn't be easier than with the flexible clamping mechanism of our RP643 series. The radiation protection for your patients is perfect and the fit is just right.

In producing our gonad aprons, we use only highquality, flexible standard lead material. The best radiation protection properties are, of course, a basic requirement.

Protective Material	Standard lead
Lead Equivalent	0.50 mm Pb
Outer Cover/Colour	ComforTex [®] HPMF, ComforTex [®] HPMF Plus (design Happy only), available in the colours:
	Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy
Sizes	S - for small children, covers 20 x 25 cm (H x W) M - for children, covers 30 x 30 cm (H x W) L - for adults, covers 37 x 40 cm (H x W) XL - for adults, covers 45 x 50 cm (H x W)

MAVIG GmbH

X-Ray Protection and Medical Suspension Systems

Stahlgruberring 5 81829 Munich Germany

PO Box 82 03 62 81803 Munich Germany

 Phone
 +49 (0) 89 / 420 96 - 0

 Fax
 +49 (0) 89 / 420 96 - 200

 e-Mail
 info@mavig.com