

EN

RP648

MAVIG



Protection during CT examinations

Reduction of radiation exposure for:

- **Thyroid gland by approx. 65 %**
- **Sternum by approx. 80 %**
- **Chest and shoulder joints by approx. 90 %**



Model: RP648

Patient protection for CT scans, e.g. of the skull
Shield for thyroid, sternum, chest, and shoulder joints

Made in Germany

CE 0302


RP648 – Patient Protection for CT Scans

RP648 – The Smart, Reusable Patient Protection for CT Scans

Patient guards for CT examinations from MAVIG provide the ideal X-ray protection. The patient can put the guards on quickly and they remain securely in the desired position.

Optimally implemented in the RP648 patient CT drape: the patient guard is easy to use and position without any further aids or closure systems. An additional thyroid protector can be folded into position if required.

- **Protective Material** Standard lead
- **Lead Equivalent** 0.35 mm Pb
- **Outer Material / Colours** ComforTex® HPMF, ComforTex® HPMF Plus (design Happy only), available in the colours:
Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy



MAVIG GmbH

X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com