

EN

RP664

MAVIG

NEW IN DESIGN "HAPPY"



Models:	RP66452S	Children Protective Apron (SMALL)
	RP66452M	Children Protective Apron (MEDIUM)
	RP66452L	Children Protective Apron (LARGE)

Made in Germany

CE 0302

RP664 – Children's Front Protection Apron

This children's apron is especially designed for universal use and therefore the optimal front protection for your little patients. The generously cut, overlapping back wings provide a balanced fit.

It should be easy to put on and remove, particularly for children's X-ray protection. The model RP664 can be easily closed at the front with a variable snap buckle closure and gives enough leeway for the exact fit. Children's X-ray protection is manufactured by MAVIG with the innovative NovaLite radiation protection material to achieve the best possible weight ratio and thus impose only the most necessary physical strain on the children.

- **Protective Material** NovaLite - radiation protective material
- **Lead Equivalent** Front protection 0.50 mm Pb / Back panels 0.25 mm Pb
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Plus (design Happy only), available in the colours:
Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy
- **Size**

SMALL	Length: 59 cm, Width: 35 cm, for children 3 to 5 years old
MEDIUM	Length: 67 cm, Width: 43 cm, for children 6 to 8 years old
LARGE	Length: 84 cm, Width: 45 cm, for children 9 to 12 years old

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com