

EN

RP679

MAVIG



Sets: **RP6791E12** Ovarian protection with stretch belt RP679 for adults (set with 3 sizes)
RP6791K12 Ovarian protection with stretch belt RP679 for children (set with 3 sizes)

Made in Germany

CE 0302

RP679 – Ovarian Protection with Stretch Belt

The patient shield RP679 protects the ovaries during X-rays of standing patients.

It is also ideally suited for children as patient protection in the appropriate size set for paediatrics and for X-rays in medical diagnostics. The stretch belt* with snap buckle closure makes it easy to put on and off and ensures convenient size adjustment.

- **X-ray protective material Standard lead with a lead equivalent of 1.00 mm Pb**
- **Comfortable size adjustment by means of the stretch belt with hygienic cover and snap buckle**
- **Quick and easy to put on**
- **Set consisting of the stretch belt* and three protectors with the following dimensions:**
 for adults: Small (Blue Admiral, 16 x 11 cm), Medium (Mandarin, 18 x 13 cm), Large (Regatta, 20 x 19 cm)
 for children: Small (Orchidee, 7 x 6 cm), Medium (Lagoon, 9 x 8 cm), Large (Curacao, 10 x 10 cm)
- **Protectors can be attached to the belt with the hook & burr fastening**

* A custom size belt is available for patients with a girth greater than 128 cm.

MAVIG GmbH
 X-Ray Protection and
 Medical Suspension Systems

Stahlgruberring 5
 81829 Munich
 Germany

PO Box 82 03 62
 81803 Munich
 Germany

Phone +49 (0) 89 / 420 96 - 0
 Fax +49 (0) 89 / 420 96 - 200
 e-Mail info@mavig.com