

EN

RP689

MAVIG



TIPP Don't be too "well-meaning" when choosing the size of your patient drape. Bear in mind that the drape also has to be moved by you and the size is directly related to the weight.

Modell: **RP689** Patient Drape (Standard lead)

Made in Germany

CE 0302

RP689 – Patient Drape

This versatile patient drape is a must in radiology practice. It provides targeted protection for your patients, combined with uncomplicated handling.

You can choose the format freely, depending on your intended use. Additional equipment with hook-and-loop straps or an individual cut according to your ideas are also possible at any time.

Effective radiation protection for your patients concerns the areas of the body adjacent to the diagnostic field; so in most cases you should rather have a stock of several, somewhat smaller drapes, which you can, of course, also use in combination.

- **Protective Material** Standard lead
- **Lead Equivalent** optionally in 0.50 mm Pb or 1.00 mm Pb
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Plus (design Happy only), available in the colours:
Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy
- **Dimensions** Our standard material widths are 60 and 90 cm, and you determine the length. Of course, you can also define a cut that deviates from the standard widths of 60 and 90 cm.
- **Edging** All edges covered with seam band

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com