



EN  
**TS100\***



MAVIG

# PORTEGRA<sup>2</sup>



<b>Modelle:</b>	<b>TS1001</b>	Versatile ceiling track width 335 mm, length 2.5 m
	<b>TS1002</b>	Versatile ceiling track width 335 mm, length 2.5 m, with cable spooler (lamps only)
	<b>TS1005</b>	Versatile ceiling track width 335 mm, length 4.0 m
	<b>TS1006</b>	Versatile ceiling track width 335 mm, length 4.0 m with cable spooler (lamps only)

Made in Germany


IEC 60601-1

## MAVIG Ceiling Tracks with or without Cable Spooler

The MAVIG ceiling track has set a new standard for system support in hospitals and radiological practices. Worldwide wherever ceiling-guided accessories, such as radiation protection shields, lamps, injectors or other equipment are utilized, MAVIG ceiling tracks are used. The unique profile structure ensures the smooth running of the carriage, while adjustable cross-struts simplify the system installation.

Medical staff and patients benefit from the flexibility in the positioning of the systems. Without physical effort, the system can be positioned and, in case of emergency, quickly removed. The carriage will always glide smoothly, even after many years of hard routine use.

- Two standard lengths: 2.5 m and 4.0 m (other lengths upon request, max. 5.0 m)
- Uniquely patented aluminum profile construction ensures greater stability
- Optional cable spooler suitable to specific track lengths
- Adjustable cross-struts, for an easy installation



**MAVIG GmbH**  
 X-Ray Protection and  
 Medical Suspension Systems

Stahlgruberring 5  
 81829 Munich  
 Germany

PO Box 82 03 62  
 81803 Munich  
 Germany

Phone +49 (0) 89 / 420 96 - 0  
 Fax +49 (0) 89 / 420 96 - 200  
 e-Mail info@mavig.com