


EN

WD504

MAVIG

Series: **WD504**Mobile X-Ray Protective Shield, 0.50 mmPb or 2.00 mmPb
Protection Lead Glass Window: 2.00 mmPb

Made in EU

WD504 – Mobile X-Ray Protective Shield

Working in radiation exposed environment requires appropriate protection. Therefore this mobile shield combines upper and lower body protection.

The WD504, available in different widths, with its wide lead glass area provides maximum protection and an unimpeded view for personnel. MAVIG offers the mobile shield WD504 in two different lead equivalent values depending on your needs.

- **MAVIG lead glass with excellent light transmission**
- **Light weight aluminum frame construction**
- **Lead equivalents** 0.50 mmPb or 2.00 mmPb (other Pb values on request)
- **Lead equivalent window** 2.00 mmPb
- **Overall widths** 700 to 1700 mm / 27.6 to 66.9 in (in steps of 100 mm / 3.9 in)
- **Overall height** 1950 mm / 76.8 in
- **Overall depth** 642 mm / 25.3 in
- **Lead glass part dimensions** Widths: 576 to 1576 mm / 22.7 to 62.0 in (in steps of 100 mm / 3.9 in, depending on overall width of mobile shield)
Height: 784 mm / 30.9 in

MAVIG GmbH

X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com