

This data sheet has been developed in order to assist with taking correct measures for tailor-made aprons. Before ordering a tailor-made apron, please double check with our customer service if not a standard size of our extensive product range would be suitable.

Please always consider the following points in order to achieve the best protection:

- Measure precisely: tape should be snug, but not restrictive
- According to R6V and DIN EN 61331-3:2016 the front of the body should be covered from the lowest point of the neck down incl. knee caps
- MAVIG recommends to even measure to the point underneath the knees for max. safety
- All of the below measures are relevant, regardless if a coat, apron or costume is requested

Please fill in this chart completely.


Apron for: Lady / Gent _____
(Name, Surname)

1 _____ cm **Total Body Height**

2 _____ cm **Breast Girth**
Measured horizontally around the body at the widest point of the chest.

3 _____ cm **Waist Girth**
Measured snug around the waist, but not restrictive.
• Ladies: at the slimmest point of the waist
• Gents: slightly below the navel

4 _____ cm **Hip-/Bottom Girth**
Measured horizontally around the body at the widest point of the hip/bottom.



Points 5 - 8 only for „Made to measure“ special orders:


5a _____ cm **Length of Apron / Coat**
Measured from the middle of the shoulder to the lower edge of the apron / coat.

5b _____ cm **Length of Costume Skirt**
Measured from waist level to the lower edge of the skirt.

6 _____ cm **Width of Back**
Measured across the back from armpit to armpit, while arms are down and close to body.

7 _____ cm **Waist Level Costume Vest**
Measured from the middle of the shoulder to the waist.
• Ladies: to the slimmest point of the waist
• Gents: to a point slightly below the navel

8 _____ cm **Circumference Upper Arm (Sleeve)**
Horizontal circumference in the middle of the upper arm.



I Model RA _____ NovaLite
 Standard Lead

Pb: Front _____ mm **Back** _____ mm
(Standard: Pb 0.25 mm)

III Sleeve right / left / both sides (optional)
 sewn-on / detachable

in Pb: _____ mm
(Standard: Pb 0.25 mm)

II Outer Material ComforTex® HPMF Hybrid
 ComforTex® HPMF

In the Colour: _____
(Colour Name or Colour Code)

IV Embroidery/Border permanently sewn-on (optional)
 detachable

Text: _____

(Please write legibly or in capital letters.)

Measured By: _____ **Date:** _____
(Name)

By/From: Hospital / Supplier / MAVIG Employee