

EN

BR115

Design Black

MAVIG



Model: BR115

Glasses with front and side X-ray protection of Pb 0.75 mm, universal size, weight approx. 110 g

Made in Germany



BR115 – X-Ray Protective Glasses

The X-ray protective glasses BR115 have slightly curved glasses and a protection of Pb 0.75 mm. The frame is made of an extremely durable plastic that ensures easy cleaning and good resistance to chemicals (disinfectant agents). The BR115 comes with a storage case and a soft google strap.

MAVIG uses only protective lenses of the highest quality. This guaranties the best light transmission values and a maximum lens break resistance. An ideal composition of the protective material results in a high level X-ray absorption, and optimized production processes achieve an absolute homogeneous allocation of the glass additives.

The quality features of our glasses make them perfect for the use in the medical environment:

- Great frame stability through high-grade synthetic materials
- Optimized balance of sturdy, light-weight frame and lenses
- Minimized risk of injury due to shape and used components
- Low allergy risk due to carefully chosen materials
- Easy cleaning and superior resistance against chemical influences (disinfection)
- Unfortunately, corrective lenses are not available for this model.
- Conformity: DIN EN 61331-1:2016 / DIN EN 61331-3:2016

MAVIG GmbH

X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com