

EN

# BR118

**Designs** Taupe  
Blue  
Black

MAVIG



Taupe



Blue



Black

**Model:** BR118

Glasses with front X-ray protection of Pb 0.75 mm, universal size, weight approx. 60 g

Made in Germany



## BR118 – X-Ray Protective Glasses

The X-ray protective glasses BR118 have curved glasses and a front protection of Pb 0.75 mm. The frame is made of an extremely durable plastic that ensures easy cleaning and good resistance to chemicals (disinfectant agents). The BR118 comes with a storage case and a soft google strap.

MAVIG uses only protective lenses of the highest quality. This guaranties the best light transmission values and a maximum lens break resistance. An ideal composition of the protective material results in a high level X-ray absorption, and optimized production processes achieve an absolute homogeneous allocation of the glass additives.

**The quality features of our glasses make them perfect for the use in the medical environment:**

- Great frame stability through high-grade synthetic materials
- Optimized balance of sturdy, light-weight frame and lenses
- Minimized risk of injury due to shape and used components
- Low allergy risk due to carefully chosen materials
- Easy cleaning and superior resistance against chemical influences (disinfection)
- For this model prescription lenses are available: single, bifocal, and progressive
- Conformity: DIN EN 61331-1:2016 / DIN EN 61331-3:2016

**MAVIG GmbH**  
X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone +49 (0) 89 / 420 96 - 0  
Fax +49 (0) 89 / 420 96 - 200  
e-Mail [info@mavig.com](mailto:info@mavig.com)