

EN

BR310

Design Tortoise

MAVIG



Model: BR310

Glasses with front and side X-ray protection of Pb 0.75 mm, universal size, weight approx. 90 g

Made in Germany



BR310 – X-Ray Protective Glasses

The X-ray protective glasses BR310 have slightly curved glasses and a protection of Pb 0.75 mm. The frame is made of an extremely durable plastic that ensures easy cleaning and good resistance to chemicals (disinfectant agents). The BR310 comes with a storage case and a soft head band.

MAVIG uses only protective lenses of the highest quality. This guaranties the best light transmission values and a maximum lens break resistance. An ideal composition of the protective material results in a high level X-ray absorption, and optimized production processes achieve an absolute homogeneous allocation of the glass additives.

The quality features of our glasses make them perfect for the use in the medical environment:

- Great frame stability through high-grade synthetic materials
- Optimized balance of sturdy, light-weight frame and lenses
- Minimized risk of injury due to shape and used components
- Low allergy risk due to carefully chosen materials
- Easy cleaning and superior resistance against chemical influences (disinfection)
- Wide and flexible head band for a comfortable fit
- For this model prescription lenses are available: single, bifocal, and progressive
- Conformity: DIN EN 61331-1:2016 / DIN EN 61331-3:2016

MAVIG GmbH

X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com