

EN

# GD4210



**Curved Support**

for one monitor, with handle



GD4210 Rear view

|                     |                  |  |
|---------------------|------------------|--|
| <b>Article-No.:</b> | <b>GD4210</b>    | Monitor support for one display up to 19", with VESA 75/100 adapter    |
|                     | <b>GD4210-L</b>  | Monitor support for one display up to 21/24", with VESA 75/100 adapter |
|                     | <b>GD4210-XL</b> | Monitor support for one display up to 27", with VESA 75/100 adapter    |

Made in Germany




## GD4210 – Single Monitor Support\* for Portegra2 Systems

The GD4210 can easily hold one monitor up to 27" (depending on version). The large curved handlebar protects the monitor and provides a means for quick & easy positioning. Combined with the mobility and safety of the Portegra2 suspension arm, the GD4210 offers the user a worry-free and flexible application.

- **Application** Monitor support with handle, VESA 75/100 adapter for one display/monitor, external cabling
- **Load / Weight**

|           |                                |
|-----------|--------------------------------|
| GD4210    | max. 15.5 kg (33 lbs) / 2.5 kg |
| GD4210-L  | max. 15.4 kg (33 lbs) / 2.6 kg |
| GD4210-XL | max. 15.2 kg (33 lbs) / 2.8 kg |
- **VESA adapter** VESA 75/100 compatible
- **Colouring** RAL9002
- **Mounting** Portegra2 systems for Stationary / ceiling track mounting
- **System range** Extension arm 750 mm / 950 mm  
Spring arm 910 mm



\* Monitor support without extension / spring arm, cabling and monitor.

\* (Dimensions for L and XL versions)

**MAVIG GmbH**  
X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone +49 (0) 89 / 420 96 - 0  
Fax +49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com