


EN

OT50001
OT90001

MAVIG



STEA-OT5 (accessory) Disposable, sterile cover for the whole shield

Article Codes: OT50001
OT90001
E-OT25B05

Shield complete with Portegra2 Suspension Arm (extension arm length 75 cm)
Shield complete with Portegra2 Suspension Arm (extension arm length 95 cm)
Replacement Shield and Connection Bar

IEC 60601-1

Made in Germany




OT50001/OT90001 – X-Ray Protective Shield

MAVIG's products are designed to offer solid protection and optimized functionality. The high quality acrylic shield allows for an excellent transmittance and unimpeded view. All edges and corners have been rounded and shatter-proof materials are used to minimize risk of injury.

MAVIG's patented Portegra2 suspension arm provides its industry known safety features and movement flexibility. The centrally guided shield with a special patient cut-out can be easily positioned over the patient and quickly removed in case of an emergency.

- **Lead Acrylic Shield** WxH: 60 x 76 cm (23.6 x 29.9 in)
- **Lead Equivalent** Pb 0.50 mm
- **OT50001** Extension / Spring Arm
Lengths: 75 cm / 91 cm (29.5 in / 35.8 in)
Colour: RAL 9002
- **OT90001** Extension / Spring Arm
Lengths: 95 cm / 91 cm (37.4 in / 35.8 in)
Colour: RAL 9002



MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com