

EN

RA614

MAVIG

Snap Fasteners

for a secure attachment of the thyroid and sternum protection to MAVIG protective clothing



Exchangeable Hook & Burr Fastener



Hygienic Cover in Colour Ocean

Model Variants: RA614KL
RA614AK
RA614MA

Thyroid/Sternum Shield with hook & burr closure
Thyroid/Sternum Shield similar to RA614KL with an exchangeable hook & burr closure
Thyroid/Sternum Shield similar to RA614KL with magnetic closure *1

Made in Germany



RA614 – Thyroid Gland and Sternum Protection

The thyroid is a radiation-sensitive organ of the human body and therefore in particular need of protection. Wearing appropriate protection during dose-intensive examinations should therefore always be a matter of course.

Our thyroid and sternum protector is characterised by its first-class cut and the best fit, taking into account all aspects of radiation protection.

With a permanently sewn-on hook & burr closure in the neck as standard, the thyroid and sternum protection is optionally also available with an exchangeable hook & burr closure or a practical magnetic closure *1.

- **Protective Material** NovaLite X-ray protective material
- **Lead Equivalent** Pb 0.50 mm
- **Outer Material / Colours** ComforTex® HPMF or ComforTex® HPMF Hybrid, available in the colours: Curacao, Limette, Orchidee, Ocean, Regatta, Lagoon, Indian Summer
- **Neck Sizes** Standard size fits most neck sizes:
RA614KL & RA614MA: 32 - 45 cm / RA614AK: 37 - 50 cm
- **Hygienic Cover** Washable hygienic cover made of the special ComforTex® MF material in the colour "Ocean" (article code: RA614HF)

(*1 not available in USA, Canada and Japan)

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com