


EN

# RD644

For Images When Sitting



Reversible Dental Apron  
available as Single Colour

or Dual Colour Design



Loops on the back to fix the RD644 to the dental chair

Article Codes: **RD644E**  
**RD644K**

X-Ray Protective Dental Apron with Collar (Thyroid Protection) for Adults  
X-Ray Protective Dental Apron with Collar (Thyroid Protection) or Children

Made in Germany



## RD644 – X-Ray Protective Dental Apron with Collar

MAVIG's products offer solid protection and optimized functionality. The dental apron model RD644 is used for images when sitting and offers high flexibility and wearing comfort. RD644 protects the whole shoulder area down to the gonads.

The practical hook & burr fastener at the thyroid protective collar ensures that it can be put on securely and quickly, with an exact fit. Loops on the back make it possible to fix the RD644 to the dental chair.

Model RD644 is designed as a reversible X-ray dental apron for patients. The design allows you to choose between a single-colour or two-colour model.

- **Protective Material** Standard lead
- **Lead Equivalent** Pb 0.50 mm
- **Outer Material / Colours** ComforTex® HPMF, available in the colours\*: Curacao, Limette, Orchidee, Ocean, Regatta, Lagoon, Indian Summer, Teddy design
- **Sizes/Dimensions** Adult: Front: 90.0 x 60.0 cm (L x W, not including collar)  
Children: Front: 70.0 x 45.0 cm (L x W, not including collar)
- **Closure** Practical hook & burr fastener at the neck
- **Loops** On the back, to fix the RD644 to the dental chair

\* Also available in a two-colour design: "Curacao with Orchidee" or "Limette with Ocean".

**MAVIG GmbH**  
X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone +49 (0) 89 / 420 96 - 0  
Fax +49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com