



RP664

MAVIG



Article No.: **RP664S**CC***** Children Protective Apron (SMALL)
RP664MCC***** Children Protective Apron (MEDIUM)
RP664LCC***** Children Protective Apron (LARGE)

Made in Germany



RP664 – Children Protective Apron

This children apron is specially designed for universal use and is thus the optimal front protection for your young patients. The generously cut, overlapping back wings offer additional protection for the back of the body and ensure a perfect fit.

Ease of putting on and taking off is particularly important in children X-ray protection. The RP664 model is easy to fasten in front with a variable snap lock and also provides enough room for exact adjustment.

- **Protective Material** NovaLite radiation protective material
- **Lead Equivalent** Front protection Pb 0.50 mm / Back panels Pb 0.25 mm
- **Outer Cover/Colour** ComforTex® HPMF, available in CC*
- **Size**

SMALL	Length: 59 cm, Width: 35 cm for children 3 to 5 years old
MEDIUM	Length: 67 cm, Width: 43 cm for children 6 to 8 years old
LARGE	Length: 84 cm, Width: 45 cm for children 9 to 12 years old

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com

(*CC - ColourCode: 1 = Curacao, 2 = Orchidee, 3 = Regatta, 4 = Indian Summer, 5 = Limette, 6 = Teddy, 7 = Ocean, 8 = Lagoon)