



RP668

MAVIG



Article No.: **RP668SCC*** Children's X-Ray Protective Coat (SMALL)
RP668MCC* Children's X-Ray Protective Coat (MEDIUM)
RP668LCC* Children's X-Ray Protective Coat (LARGE)

Made in Germany



RP668 – Children's X-Ray Protective Coat

With this model, you have chosen the practical all around protection for your little patients. The cut combines ideal radiation protection with a variable fit and a high degree of wearing comfort. At the same time the coat is easy to put on and take off.

The RP668 comes with hook & burr fasteners on the front for a quick and safe adjustment of the coat. MAVIG produces children's X-ray protection with the innovative NovaLite radiation protection material in order to achieve the best possible weight ratio. Thus, children only have to endure the lowest possible weight.

- **Protective Material** NovaLite radiation protective material
- **Lead Equivalent** Front protection Pb 0.50 mm / Back panels Pb 0.25 mm
- **Outer Cover/Colour** ComforTex® HPMF, available in CC*
- **Size**

SMALL	Length: 60 cm, Hip girth: 68 - 75 cm for children 3 to 5 years old
MEDIUM	Length: 80 cm, Hip girth: 76 - 83 cm for children 6 to 8 years old
LARGE	Length: 84 cm, Hip girth: 84 - 95 cm for children 9 to 12 years old

MAVIG GmbH
 X-Ray Protection and
 Medical Suspension Systems

Stahlgruberring 5
 81829 Munich
 Germany

PO Box 82 03 62
 81803 Munich
 Germany

Phone +49 (0) 89 / 420 96 - 0
 Fax +49 (0) 89 / 420 96 - 200
 e-Mail info@mavig.com

(*CC - ColourCode: 1 = Curacao, 2 = Orchidee, 3 = Regatta, 4 = Indian Summer, 5 = Limette, 6 = Teddy, 7 = Ocean, 8 = Lagoon)