


EN RP270

MAVIG



Article No.: RP270 Breast Protection for Women

Made in Germany




RP270 – Breast Protection for Women

MAVIG's breast protection was specially designed for the protection of the mammary gland area of the patients during mammographic radiation procedures. As all of our protection the RP270 offers solid protection and optimized functionality.

The protector consists of two parts, each 18 x 40 cm (W x H), joined by an adaptable hook & burr closure for a comfortable fit and easy size adjustment.

- **Protective Material** Novalite radiation protective material
- **Lead Equivalent** Pb 0.50 mm
- **Size** Universal
- **Outer Cover/Colour** ComforTex® HPMF, available in the colours: *Curacao, Limette, Regatta, Indian Summer, Ocean, Lagoon, Orchidee and Teddy Design*
- **Edging** All edges are covered with a seem band.



MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com