

With the help of this data sheet and a standard tape measure, you can easily determine your body measurements. Based on the measurements, a clothing size will be determined for you or, if necessary, a made-to-measure item will be suggested. Before you use the measuring tape, please note the following:

The table data are to be understood as body measurements, measured **without additions** but also without being constricting. You **cannot** determine the correct values if you measure yourself. The measurements must thus always be taken by another person. The person being measured should stand straight with arms outstretched.

The length of the X-ray protection clothing should be selected in accordance with the German Radiation Protection Ordinance (RöV) and DIN EN 61331-3:2016 so that the front of the body is covered from the base of the neck to at least the kneecaps (MAVIG recommends coverage to below the knee for complete protection).

**Please fill in the following points completely.**


**Measurements for:**  Mr /  Ms \_\_\_\_\_  
(First name and surname)

**1** \_\_\_\_\_ cm **Total height**

**2** \_\_\_\_\_ cm **Chest circumference**  
Across the broadest part of the chest (bust) horizontally around the body

**3** \_\_\_\_\_ cm **Waist circumference**  
Measure around the waist without making it tight (for women measure at the narrowest part of the waist, for men measure around the body slightly below the navel)

**4** \_\_\_\_\_ cm **Hip/buttock circumference**  
Place the measuring tape horizontally around the broadest part of the hips/buttocks



**Points 5 – 8 only for custom-made products:**


**5a** \_\_\_\_\_ cm **Apron/coat length**  
Measure from the middle of the shoulder, down the chest to the desired height of the lower edge of the apron/coat

**5b** \_\_\_\_\_ cm **Costume skirt length**  
Measure from the waist down to the desired height of the bottom edge of the skirt

**6** \_\_\_\_\_ cm **Back width**  
Measure at the back, from armpit to armpit, with arms at the sides

**7** \_\_\_\_\_ cm **Front length of the costume vest**  
Measure from the middle of the shoulder over the chest down to the narrowest part of the waist (for women) or to slightly below the navel (for men)

**8** \_\_\_\_\_ cm **Upper arm circumference (sleeve)**  
Horizontal circumference at the middle of the upper arm



**I Model** RA  NovaLite  
 Leadfree  
 Standard lead

**Pb: Front** \_\_\_\_\_ mm **Back** \_\_\_\_\_ mm  
(Standard: 0.25 mmPb)

**III Sleeves**  Right /  Left /  Both sides  
(Optional)  Sewn-on /  Buttoned-on

**Pb:** \_\_\_\_\_ mm  
(Standard: 0.25 mmPb)

**II Outer material**  ComforTex® HPMF Hybrid  
 ComforTex® HPMF  
 ComforTex® HPMF Plus  
(Nautilus design only)

**Colour:** \_\_\_\_\_  
(Colour indication or colour code)

**IV Embroidery/border**  Sewn-in  
(Optional)  Adjustable

**Text:** \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
(Please write legibly, in capital letters)

**Measured by:** \_\_\_\_\_ **on:** \_\_\_\_\_  
(Name) (Date)

**From/by:**  Hospital /  Supplier /  MAVIG employee