



Article-No.:	AW42512	Mobile storage system
	AW425K	Basket for accessories, compatible with AW42512

AW425 – Mobile Storage System for X-Ray Protective Clothing

Storage systems for X-ray protective clothing must fulfil a number of basic criteria, in order to be of reliable use in your daily work.

With the AW425, your X-ray protective clothing is safely stored using the MAVIG special hangers, included with every apron free of charge. The smooth-running storage system offers storage space for up to twelve aprons.

- **Convenient & simple access**
- **Clear arrangement for up to 12 pieces of clothing**
- **Minimal space requirement**
- **Suitable and gentle storage**
- **Easy assembly**

• Technical Data	Weight (net)	15.2 kg / basket 2.0 kg
	Mobility	4 heavy load casters, 2 with brakes
	Storage capacity	12 special hangers / X-ray protective aprons (max. load capacity 130.0 kg)
	Outer dimensions	1000 x 1540 x 600 mm (W x H x D)
	Colour	Stainless steel / RAL 9010

MAVIG GmbH

X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com