

EU-Konformitätserklärung *EU Declaration of Conformity*

Name und Anschrift des Herstellers/
Name and Address of Manufacturer: MAVIG GmbH
Stahlgruberring 5
81829 Munich
Germany

Für die Ausstellung dieser EU-Konformitätserklärung trägt MAVIG GmbH die alleinige Verantwortung und bestätigt, dass die nachfolgend genannten Produkte die Anforderungen der Verordnung (EU) 2016/425 für Persönliche Schutzausrüstung erfüllen.

This EU declaration of conformity is issued under the sole responsibility of MAVIG GmbH who confirms that the products mentioned below meet the provision of the Regulation (EU) 2016/425 on personal protective equipment.

Produktgruppe/
Product Group: RA Röntgenschutzkleidung Anwender
RA Radiation Protective Clothing Medical

Produktbezeichnung/
Product name: Rundumschutz-Mantel Victoria
Coat Victoria

Modell/ Model: RA636-STD (Standard/*standard*)
RA636-SON (Sonderanfertigung/*custom made*)
RA636-LVO (Levio Edition)

Angewandte Normen/
Applied standards: IEC 61331-1:2014
IEC 61331-3:2014

Die benannte Stelle ANCCP Certification Agency srl., Via dello Struggino, 6 – 57121 Livorno (LI) – Italy hat die EU-Baumusterprüfung durchgeführt und die EU-Baumusterprüfbescheinigung PPE-2094-21811-Cert (PPE-2095-21811-Cert für Levio Edition) ausgestellt.

The notified body ANCCP Certification Agency srl., Via dello Struggino, 6 – 57121 Livorno (LI) – Italy performed the EU type-examination and issued the EU type-examination certificate PPE-2094-21811-Cert (PPE-2095-21811-Cert for Levio Edition).

Die PSA unterliegt dem in Anhang VIII (Modul D) der Verordnung (EU) 2016/425 vorgesehenen Verfahren der regelmäßigen Überwachung unter der Aufsicht der folgenden Benannten Stelle:

The PPE is subject to the scheduled periodic surveillance procedure set out in Annex VIII (Module D) of Regulation (EU) 2016/425 under the supervision of the following Notified Body:

ANCCP Certification Agency srl., Via dello Struggino, 6 – 57121 Livorno (LI) – Italy – CE0302

Diese Erklärung wurde ausgestellt von/
This statement was issued by:

Ort und Datum der Ausstellung/
Place and date of issue:

Munich, 2024.01.30

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Company Address:
Stahlgruberring 5
81829 Munich
Germany

Post Box:
Postfach 820362
81803 Munich
Germany

Phone:
+49 (0) 89 / 420 96 - 0

Fax:
+49 (0) 89 / 420 96 - 200

e-Mail:
info@mavig.com

Web:
www.mavig.com

Bank Accounts:
HypoVereinsbank Munich
IBAN: DE80 7002 0270
1620 2253 00
SWIFT-BIC: HYVEDEMMXXX

Postbank Munich
IBAN: DE12 7001 0080
0071 3748 07
SWIFT-BIC: PBNKDEFFXXX

Trade-Register Munich
HRB-No. 76 403

Managing Directors:
Christian Stoian
Alexia Stoian-Lepère
Philip Milbredt

VAT-ID No.:
DE 811268845

Modellvarianten/ *Model versions:*

Beispielhafte Zusammensetzung des Artikelcodes/
Example for the formation of the article code:

RA636	3	2	Z	65	22	2	8
	①	②	③	④	⑤	⑥	⑦

① Bleigleichwert-Front/ Lead Equivalent-Front	
	Code
0.25 mm Pb	2
0.35 mm Pb	3
0.50 mm Pb	5
0.25 mm Pb / 0.42 mm Pb (Levio)	6
0.25 mm Pb / 0.50 mm Pb (Levio)	7

② Bleigleichwert-Rücken/ Lead Equivalent-Back	
	Code
0.25 mm Pb	2

③ Größen/ Sizes			
	Code		Code
XS	Z	XL	X
S	S	XXL	D
M	M	3XL	T
L	L		

④ Längen/ Lengths	
	Code
100 cm	10
110 cm	11
120 cm	12
130 cm	13
140 cm	14
90 cm	90

⑤ Farbauswahl/ Colours			
	Code		Code
Curacao	22	Titan	29
Lagoon	25	Mandarin	30
Orchidee	27	Blue Admiral	31
Regatta	28	Baccara	32

⑥ Außenmaterial – Hülle/ Outer Material - Cover	
	Code
ComforTex HPMF	1
ComforTex HPMF Hybrid	2

⑦ Strahlenschutzmaterial – Innen/ X-Ray Protective Material - Inlay	
	Code
NovaLite 0.25 mm Pb and 0.50 mm Pb	1
NovaLite 0.35 mm Pb	8
Standard-Blei/ Standard Lead 0.25 mm Pb, 0.35 mm Pb and 0.50 mm Pb	2
Leadfree	3

MAVIG GmbH
 X-Ray Protection and
 Medical Suspension Systems

Trade-Register Munich
 HRB-No. 76 403

Managing Directors:
 Christian Stoian
 Alexia Stoian-Lepère
 Philip Milbrdt

VAT-ID No. :
 DE 811268845