

Modellvarianten/ Model versions:

Beispielhafte Zusammensetzung des Artikelcodes/
Example for the formation of the article code:

RA631	K	3	2	Z	65	22	2	8
	①	②	③	④	⑤	⑥	⑦	⑧

① Ausführung/ Type	
	Code
Kostüm/ Costume	K
Weste/ Vest	W
Rock/ Skirt	R
Weste Victoria/ Vest Victoria	V

② Bleigleichwert-Front/ Lead Equivalent-Front	
	Code
0.25 mm Pb	2
0.35 mm Pb	3
0.50 mm Pb	5
0.25 mm Pb / 0.42 mm Pb (Levio)	6
0.25 mm Pb / 0.50 mm Pb (Levio)	7

③ Bleigleichwert-Rücken/ Lead Equivalent-Back	
	Code
0.25 mm Pb	2

④ Größen/ Sizes			
	Code		Code
XS	Z	XL	X
S	S	XXL	D
M	M	3XL	T
L	L		

⑤ Rocklängen/ Skirt Length	
	Code
65 cm	65
70 cm	70
75 cm	75
80 cm	80

⑤ Farbauswahl/ Colours			
	Code		Code
Curacao	22	Titan	29
Lagoon	25	Mandarin	30
Orchidee	27	Blue Admiral	31
Regatta	28	Baccara	32

⑥ Außenmaterial – Hülle/ Outer Material - Cover	
	Code
ComforTex HPMF	1
ComforTex HPMF Hybrid	2

⑦ Strahlenschutzmaterial – Innen/ X-Ray Protective Material - Inlay	
	Code
NovaLite 0.25 mm Pb and 0.50 mm Pb	1
NovaLite 0.35 mm Pb	8
Standard-Blei/ Standard Lead 0.25 mm Pb, 0.35 mm Pb and 0.50 mm Pb	2
Leadfree	3

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Trade-Register Munich
HRB-No. 76 403

Managing Directors:
Christian Stoian
Alexia Stoian-Lepère
Philip Milbrdt

VAT-ID No.:
DE 811268845