

## EU-Konformitätserklärung *EU Declaration of Conformity*

Name und Anschrift des Herstellers/  
*Name and Address of Manufacturer:*

MAVIG GmbH  
Stahlgruberring 5  
81829 Munich  
Germany

Für die Ausstellung dieser EU-Konformitätserklärung trägt MAVIG GmbH die alleinige Verantwortung und bestätigt, dass die nachfolgend genannten Produkte die Anforderungen der Verordnung (EU) 2016/425 für Persönliche Schutzausrüstung erfüllen.

*This EU declaration of conformity is issued under the sole responsibility of MAVIG GmbH who confirms that the products mentioned below meet the provision of the Regulation (EU) 2016/425 on personal protective equipment.*

Produktgruppe/  
*Product Group:*

RA Röntgenschutzkleidung Anwender  
*RA Radiation Protective Clothing Medical*

Produktbezeichnung/  
*Product name:*

Frontschuttschürze Basic  
*Front Apron Basic*

Modell/ *Model:*

RA640-STD (Standard/*standard*)  
RA640-SON (Sonderanfertigung/*custom made*)  
RA640-LVO (Levio Edition)

Angewandte Normen/  
*Applied standards:*

IEC 61331-1:2014  
IEC 61331-3:2014

Die benannte Stelle ANCCP Certification Agency srl., Via dello Struggino, 6 – 57121 Livorno (LI) – Italy hat die EU-Baumusterprüfung durchgeführt und die EU-Baumusterprüfbescheinigung PPE-2094-21811-Cert (PPE-2095-21811-Cert für Levio Edition) ausgestellt.

*The notified body ANCCP Certification Agency srl., Via dello Struggino, 6 – 57121 Livorno (LI) – Italy performed the EU type-examination and issued the EU type-examination certificate PPE-2094-21811-Cert (PPE-2095-21811-Cert for Levio Edition).*

Die PSA unterliegt dem in Anhang VIII (Modul D) der Verordnung (EU) 2016/425 vorgesehenen Verfahren der regelmäßigen Überwachung unter der Aufsicht der folgenden Benannten Stelle:

*The PPE is subject to the scheduled periodic surveillance procedure set out in Annex VIII (Module D) of Regulation (EU) 2016/425 under the supervision of the following Notified Body:*

ANCCP Certification Agency srl., Via dello Struggino, 6 – 57121 Livorno (LI) – Italy – CE0302

Diese Erklärung wurde ausgestellt von/  
*This statement was issued by:*



Ort und Datum der Ausstellung/  
*Place and date of issue:*

Munich, 2024.01.30

MAVIG GmbH  
X-Ray Protection and  
Medical Suspension Systems

Company Address:  
Stahlgruberring 5  
81829 Munich  
Germany

Post Box:  
Postfach 820362  
81803 Munich  
Germany

Phone:  
+49 (0) 89 / 420 96 - 0

Fax:  
+49 (0) 89 / 420 96 - 200

e-Mail:  
info@mavig.com

Web:  
www.mavig.com

Bank Accounts:  
HypoVereinsbank Munich  
IBAN: DE80 7002 0270  
1620 2253 00  
SWIFT-BIC: HYVEDEMMXXX

Postbank Munich  
IBAN: DE12 7001 0080  
0071 3748 07  
SWIFT-BIC: PBNKDEFFXXX

Trade-Register Munich  
HRB-No. 76 403

Managing Directors:  
Christian Stoian  
Alexia Stoian-Lepère  
Philip Milbredt

VAT-ID No.:  
DE 811268845

Modellvarianten/ *Model versions:*

Beispielhafte Zusammensetzung des Artikelcodes/  
*Example for the formation of the article code:*

|              |   |   |   |    |    |    |   |   |
|--------------|---|---|---|----|----|----|---|---|
| <b>RA640</b> | 3 | 2 | Z | 12 | 22 | DK | 2 | 8 |
|              | ① | ② | ③ | ④  | ⑤  | ⑥  | ⑦ | ⑧ |

| ①<br>Bleigleichwert-Front/<br>Lead Equivalent-Front |      |
|---|------|
|   | Code |
| 0.25 mm Pb  | 2    |
| 0.35 mm Pb  | 3    |
| 0.50 mm Pb  | 5    |
| 0.25 mm / 0.42 mm Pb (Levio)                        | 6    |
| 0.25 mm / 0.50 mm Pb (Levio)                        | 7    |

| ②<br>Bleigleichwert-Rücken/<br>Lead Equivalent-Back |      |
|---|------|
|   | Code |
| 0.25 mm Pb  | 2    |
| Ohne Schutz/ Without protection                     | 0    |

| ③<br>Größen/ Sizes |      |     |      |
|--------------------|------|-----|------|
|                    | Code |     | Code |
| XS                 | Z    | XL  | X    |
| S                  | S    | XXL | D    |
| M                  | M    | 3XL | T    |
| L                  | L    |     |      |

| ④<br>Längen/ Lengths |      |
|----------------------|------|
|                      | Code |
| 100 cm               | 10   |
| 110 cm               | 11   |
| 120 cm               | 12   |
| 130 cm               | 13   |
| 140 cm               | 14   |
| 90 cm                | 90   |

| ⑤<br>Farbauswahl/ Colours |      |              |      |
|---------------------------|------|--------------|------|
|                           | Code |              | Code |
| Curacao                   | 22   | Titan        | 29   |
| Lagoon                    | 25   | Mandarin     | 30   |
| Orchidee                  | 27   | Blue Admiral | 31   |
| Regatta                   | 28   | Baccara      | 32   |

| ⑥<br>Verschlussart/ Type of Closure |      |
|-------------------------------------|------|
|                                     | Code |
| Druckknopf/ Snap Fasteners          | DK   |
| Kreuzband/ Crossbands               | KB   |
| Gurtband/ Straps & Snap Buckle      | GB   |
| Schnappschnallen/ Snap Buckles      | SV   |

| ⑦<br>Außenmaterial – Hülle/<br>Outer Material - Cover |      |
|---|------|
|   | Code |
| ComforTex HPMF  | 1    |
| ComforTex HPMF Hybrid                                 | 2    |

| ⑧<br>Strahlenschutzmaterial – Innen/<br>X-Ray Protective Material - Inlay |      |
|---|------|
|   | Code |
| Novalite 0.25 mm Pb and 0.50 mm Pb  | 1    |
| Novalite 0.35 mm Pb   | 8    |
| Standard-Blei/ Standard Lead 0.25 mm Pb, 0.35 mm Pb and 0.50 mm Pb        | 2    |
| Leadfree  | 3    |

MAVIG GmbH  
 X-Ray Protection and  
 Medical Suspension Systems

Trade-Register Munich  
 HRB-No. 76 403

Managing Directors:  
 Christian Stoian  
 Alexia Stoian-Lepère  
 Philip Milibredt

VAT-ID No.:  
 DE 811268845