

Modellvarianten/ *Model versions:*

Beispielhafte Zusammensetzung des Artikelcodes/
Example for the formation of the article code:

RA641	3	2	L	12	22	GB	2	8
	①	②	③	④	⑤	⑥	⑦	⑧

① Bleigleichwert-Front/ Lead Equivalent-Front	
	Code
0.25 mm Pb	2
0.35 mm Pb	3
0.50 mm Pb	5
0.25 mm / 0.42 mm Pb (Levio)	6
0.25 mm / 0.50 mm Pb (Levio)	7

② Bleigleichwert-Rücken/ Lead Equivalent-Back	
	Code
0.25 mm Pb	2
Ohne Schutz/ Without protection	0

③ Größen/ Sizes			
	Code		Code
XS	Z	XL	X
S	S	XXL	D
M	M	3XL	T
L	L		

④ Längen/ Lengths	
	Code
100 cm	10
110 cm	11
120 cm	12
130 cm	13
140 cm	14
90 cm	90

⑤ Farbauswahl/ Colours			
	Code		Code
Curacao	22	Titan	29
Lagoon	25	Mandarin	30
Orchidee	27	Blue Admiral	31
Regatta	28	Baccara	32

⑥ Verschlussart/ Type of Closure	
	Code
Gurtband/ Straps & Snap Buckle	GB
Schnappschnallen/ Snap Buckles	SV

⑦ Außenmaterial – Hülle/ Outer Material - Cover	
	Code
ComforTex HPMF	1
ComforTex HPMF Hybrid	2

⑧ Strahlenschutzmaterial – Innen/ X-Ray Protective Material - Inlay	
	Code
Novalite 0.25 mm Pb and 0.50 mm Pb	1
Novalite 0.35 mm Pb	8
Standard-Blei/ Standard Lead 0.25 mm Pb, 0.35 mm Pb and 0.50 mm Pb	2
Leadfree	3

MAVIG GmbH
 X-Ray Protection and
 Medical Suspension Systems

Trade-Register Munich
 HRB-No. 76 403

Managing Directors:
 Christian Stoian
 Alexia Stoian-Lepère
 Philip Milbredt

VAT-ID No. :
 DE 811268845