




Standard cover version, at lowest extension arm

NEW




Hybrid Central Axis Flat Cover

Flat Cover Version – Advantages

- **Save space: 29.5 mm less in height**
(compared to standard version / height: 54.5 mm)
- **Higher flexibility for room planning**
- **Avoid conflict with other equipment**
(e.g. stationary/mobile suspension arms, movable C-arm ...)
- **Easy exchange**

Example



Article No.: HZA-COV

Complete set with flat cover (height: 25 mm) and four screws for mounting
(the set is optional, not included as standard, and only available upon request)

Beware: When using the flat cover, internal cabling on the lowest arm is not possible.



MAVIG GmbH

Stahlgruberring 5
81829 Munich
Germany

Phone +49 (0) 89 420 96 0
Fax +49 (0) 89 420 96 200
e-Mail info@mavig.com

www.mavig.com