

Tailor Made Orders

With the help of this template and a standard measuring tape, you can easily determine your body measurements.

Based on the measurements, a made-to-measure apron will be manufactured. Before you use the measuring tape, please note the following:

The requested measurements should be body measurements, measured **without any additions** but also without being constricting. Correct values cannot be determined if you measure yourself. The measurements must always be taken by another person. The person being measured should stand straight with arms hanging to the sides.

The length of the X-ray protection clothing should be selected in accordance with the German Radiation Protection Ordinance (R6V) and DIN EN 61331-3:2016 so that the front of the body is covered from the base of the neck to at least the kneecaps (MAVIG recommends coverage to below the knee for complete protection).

Please fill in the following points completely.

● **Measurements for:** _____
(First name and surname)

1 _____ cm **Total height**

2 _____ cm **Chest circumference**

Across the broadest part of the chest (bust) horizontally around the body

3 _____ cm **Waist circumference**

Measure around the waist without making it tight (for women measure at the narrowest part of the waist, for men measure around the body slightly below the navel)

4 _____ cm **Hip/buttock circumference**

Place the measuring tape horizontally around the broadest part of the hips/buttocks

5a _____ cm **Apron/coat length**

Measure from the middle of the shoulder, down the chest to the desired height of the lower edge of the apron/coat

5b _____ cm **Costume skirt length**

Measure from the waist down to the desired height of the bottom edge of the skirt

6 _____ cm **Back width**

Measure at the back, from armpit to armpit, with arms at the sides

7 _____ cm **Front length of the upper body**

Measure from the middle of the shoulder over the chest down to the Iliac crest (Spinae iliacae anteriores)

8 _____ cm **Upper arm circumference (sleeves, optional)**

Horizontal circumference at the middle of the upper arm

9 _____ cm **Position from above (only Levio®)**



Measure from the middle of the shoulder the start of the breast (about at the second rib, between the collarbone and the mammary papilla)

10 _____ cm **Length / Skirt panel length (only Levio®)**

a) Measure from the middle of the shoulder to below the gonads (about 10 cm below the crotch)

b) Skirt only: Measure from the waist (waistband of the skirt) to below the gonads (about 10 cm below the crotch)

Woman Man



I **Model** RA _____ Levio Edition
 Victoria

Inside: NovaLite Plus Leadfree 2.0

Pb Front: _____ mm **Back:** _____ mm
(Standard: 0.25 mm Pb)

III **Sleeves** Right / Left / Pair
(optional) Sewn-on / Buttoned-on

in Pb: _____ mm
(Standard: 0.25 mm Pb)

Bolero in Pb: _____ mm
(Standard: 0.25 mm Pb)

II **Outer Material** ComforTex® HPMF Hybrid
 ComforTex® HPMF

Optional: with Outlast®, optional

Outer Colour: _____
(Colour name or colour code)

IV **Embroidery/Border** Sewn-in Adjustable
(optional) Shoulder Skirt

Text: _____

(Please write legibly, in capital letters)

● **Measured by:** _____ **on:** _____
(Name) (Date)

From/by: Hospital / Supplier / MAVIG employee