

EN

RA637 Bolero

MAVIG



Model: RA637 Shoulder protector "Bolero", matching all MAVIG all-round protection or front protection aprons

Made in Germany

RA637: Shoulder Protector – „Bolero“

Thanks to the newly developed Bolero as an alternative to the sleeve, the radiation exposure due to scattered radiation for the shoulder joint can be significantly reduced. Freedom of movement is maintained to the greatest extent possible with the design, which is open underneath the arm. At the same time, the design considerably reduces the scattered radiation coming in from the sides via the armholes.

The "Bolero" shoulder protector is available for all models and is therefore ideal for combination with MAVIG all-round protection or front protection aprons. It is worn directly under the X-ray protection and is simply closed with a hook-and-loop fastener at the front.

The standard lead value is 0.25 mm Pb. This ensures good mobility during the daily working routine. The shoulder protector is always manufactured in the colour "Titan". The correct size as well as the protective material depend on the chosen X-ray protective clothing.

- **Protects the AC shoulder joint with standard: 0.25 mm Pb**
(on request: 0.35 mm Pb or 0.50 mm Pb)
- **worn directly under the X-ray protection**
- **Protective material options: NovaLite® Plus or Leadfree 2.0**
- **Outer material in the colour Titan, optional: ComforTex® HPMF or ComforTex® HPMF Hybrid**
- **Sizes ranging from XS to 3XL** (matching the PPE)

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com