

Ergonomic RA650

RA650 radiation protection apron – Comprehensive protection with maximum comfort

The RA650 offers you extensive protection for the entire front part of your body covering 60% of your body circumference. This apron is ideal for situations in which you are only exposed to radiation from the front.

Your benefits at a glance:

- ✔ **Easy Handling:** wide cross straps make it easy to remove the apron. The hook and loop fasteners at the front make donning and doffing quick and uncomplicated.
- ✔ **Lightweight:** Choose the option without protection in the back area to reduce the weight of the apron and enjoy even more wearing comfort.
- ✔ **Added Comfort:** shoulder pads are included as a standard for added comfort during wearing.
- ✔ **Long-lasting Closure Systems:** We only use high-quality hook and loop fasteners that are particularly resistant to wear and tear. These ensure a high degree of durability and easy handling of the apron.
- ✔ **Certified Protection:** Type testing and certification in accordance with (EU) 2016/425, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016



FRONT PROTECTION

Protection/Lead Equivalent Standard Collection

Novalite® Plus or Leadfree 2.0



Front:

0.25 mm Pb, 0.35 mm Pb
or 0.50 mm Pb

Back:

0.25 mm Pb (standard)
or 0.00 mm Pb (optional)



Sizes*, recommended for:

XS	Chest girth ≤ 84 cm	Waist/hip girth ≤ 91 cm
S	Chest girth approx. 85 - 95 cm	Waist/hip girth approx. 92 - 102 cm
M	Chest girth approx. 96 - 106 cm	Waist/hip girth approx. 103 - 112 cm
L	Chest girth approx. 107 - 117 cm	Waist/hip girth approx. 113 - 123 cm
XL	Chest girth approx. 118 - 128 cm	Waist/hip girth approx. 124 - 134 cm
XXL	Chest girth approx. 129 - 139 cm	Waist/hip girth approx. 135 - 145 cm
3XL	Chest girth approx. 140 - 150 cm	Waist/hip girth approx. 146 - 156 cm

* The size depends on the girth of the chest, waist, or hip, depending on which of the three measurements is the largest.

Choose the RA650 - for comprehensive front protection, high wearing comfort, and reliable quality.



Outlast Temperature Regulation



Application Oriented

The apron can be slid off easily and quickly by simply opening the hook and loop fasteners at the front. This well thought-out function allows the apron to be removed quickly and easily with minimal effort, even when worn under a sterile surgical gown.

For an application oriented use this contributes significantly to increasing efficiency and optimizing the workflow without compromising comfort or sterility.

Fast, Simple and Convenient

The apron can be slid off by simply opening the hook-and-loop fasteners in the side area, even through the sterile surgical gown.

Shoulder Relief/Crossbands

Thanks to the extra-wide, comfortable crossbands, the load is optimally distributed over the entire shoulder area.


Shoulder Pads

Also our front aprons come with soft shoulder pads that take the strain off your shoulders.



Outer Material

Colours
(see front page)



Optional:
ComforTex® HPMF,
ComforTex® HPMF Hybrid

Sizes

from XS to 3XL
(see front page)



Lengths / Body height:
90 cm / 150 - 160 cm
100 cm / 160 - 170 cm
110 cm / 170 - 180 cm
120 cm / 175 - 190 cm
130 cm / 185 - 195 cm
140 cm / 195 - 205 cm

Optional



On request:
Outlast
Breast pocket
Embroidery
Made-to-measure protection

Accessories

Bolero
Thyroid/sternum protection

