

# Comfort RA660

## RA660 front protection apron – Modern protection for your everyday work

Our customers value the RA660 as their favorite among the front protection aprons. This apron is specially designed for applications where you are only exposed to radiation from the front. It is easy to put on and take off, adapts nicely to the body and offers an excellent fit and shoulder relief.

Why the RA660 is the best choice for you:

- ✔ **Outstanding Comfort:** the wide stretch inserts, which extend over the entire back, ensure that the apron fits well and provides an optimum fit, especially when worn for long periods of time.
- ✔ **Flexibility:** Choose the optional version without protection in the back area to reduce the weight of the apron and enjoy even more comfort.
- ✔ **Shoulder Relief:** The apron comes as a standard with soft shoulder pads to relieve the strain on your shoulders.
- ✔ **Closure System:** For maximum durability, we only use high-quality hook and loop fasteners that are particularly durable.
- ✔ **Certified Protection:** Type testing and certification in accordance with (EU) 2016/425, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016

### Protection/Lead Equivalent Standard Collection

NovaLite® Plus or Leadfree 2.0



**Front:**  
0.25 mm Pb, 0.35 mm Pb  
or 0.50 mm Pb

**Back generally:**  
0.25 mm Pb



### Protection/Lead Equivalent Levio Edition

NovaLite® Plus or Leadfree 2.0



**Total value Levio-Panel:**  
0.42 mm Pb or 0.50 mm Pb  
(outside 0.25 mm Pb)

**Back generally:**  
0.25 mm Pb



Sizes\*, recommended for:

XS	Chest girth ≤ 84 cm	Waist/hip girth ≤ 91 cm
S	Chest girth approx. 85 - 95 cm	Waist/hip girth approx. 92 - 102 cm
M	Chest girth approx. 96 - 106 cm	Waist/hip girth approx. 103 - 112 cm
L	Chest girth approx. 107 - 117 cm	Waist/hip girth approx. 113 - 123 cm
XL	Chest girth approx. 118 - 128 cm	Waist/hip girth approx. 124 - 134 cm
XXL	Chest girth approx. 129 - 139 cm	Waist/hip girth approx. 135 - 145 cm
3XL	Chest girth approx. 140 - 150 cm	Waist/hip girth approx. 146 - 156 cm

\* The size depends on the girth of the chest, waist, or hip, depending on which of the three measurements is the largest.



FRONT PROTECTION

*Decide on the RA660 - for modern protection, high comfort, and outstanding quality in your day-to-day work.*



## Comfortable Design

### Wide Stretch Insert

The insert does not only ensure a secure and safe fit of the apron, but also a perfect adaptation to individual body shapes. The sophisticated design enables an optimum fit that distributes pressure evenly and noticeably relieves pressure on the shoulders and back.

This innovative solution significantly increases comfort and prevents fatigue, especially when worn for long periods of time, so that you can fully concentrate on your work.

### Levio® Edition

The apron RA660 is also available in our Levio Edition. Optimized weight-to-protection ratio: Instead of a uniform lead equivalent across the entire front area, Levio protects - from above the chest to below the gonads - with a higher lead equivalent.

### Outlast® – Ideal Temperature for Well-being

Optionally, the the apron RA660 can be equipped with Outlast technology. The Outlast® inlay material absorbs heat specifically, stores and releases it again when needed.



Outlast Temperature Regulation

### Higher Protection – Levio Edition

The area with the reinforced protective material is clearly recognisable by its colour.

### User Friendly



The sophisticated design of the RA660 front protection apron makes it easy and quick to put on or take off.

### Special Feature Stretch insert

The stretch insert offers ideal relief especially for long periods of wear. The stretch insert running across the entire width of the back ensures a good fit and shoulder relief.

### Shoulder Pads

Also our front aprons come with soft shoulder pads that take the strain off your shoulders.



### Outer Material

Colours  
(see front page)



Optional:  
Comfortex® HPMF,  
Comfortex® HPMF Hybrid

### Sizes

from XS to 3XL  
(see front page)



Lengths / Body height:  
90 cm / 150 - 160 cm  
100 cm / 160 - 170 cm  
110 cm / 170 - 180 cm  
120 cm / 175 - 190 cm  
130 cm / 185 - 195 cm  
140 cm / 195 - 205 cm

### Optional



On request:  
Outlast  
Sleeves  
Breast pocket  
Embroidery  
Made-to-measure protection

### Accessories

Bolero  
Thyroid/sternum protection

