

EN
RD614


MAVIG




Note: In order to meet hygiene requirements, we recommend that you use our washable hygienic cover made from special ComforTex® MF material "Titan" colour.



With interchangeable hook-and-loop fastener



Children's size in Happy design



Available optional hygienic cover

Models: **RD6145E** Adult patient thyroid and sternum protection
RD6145K Children patient thyroid and sternum protection

Made in Germany

CE 0302

RD614 – Thyroid and Sternum Protection for Patients

The thyroid is a radiation-sensitive organ of the human body and therefore in particular need of protection. The MAVIG thyroid sternum protection for patients is optimally cut and has the best fit, and, of course, it also takes all radiation protection aspects into account.

An uncomplicated magnetic or hook-and-loop fastener, which can also be replaced as an option, provides sufficient scope for adjustment.

- **Protective Material** Novalite Plus or Leadfree 2.0
- **Lead Equivalent** 0.50 mm Pb
- **Magnetic/hook-and-loop fastener**
 - Adults Standard size fits neck sizes of 32 - 45 cm
 - Children Standard size fits neck sizes of 28 - 36 cm
- **Interchangeable hook-and-loop fastener**
 - Adults Standard size fits neck sizes of 37 - 50 cm
 - Children Standard size fits neck sizes of 28 - 36 cm
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Hybrid, available in the colours:
Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy (ComforTex® HPMF only)

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com