

EN

RD635

MAVIG



Models: **RD6355E**
RD6355K

Adult patient X-ray apron for panorama shots while standing
Children patient X-ray apron for panorama shots while standing

Made in Germany

CE 0302

RD635 – Patient X-Ray Apron for Panorama Shots

The RD635 protective apron for patients offers optimum radiation protection for standing X-ray exposures on panoramic or teleradiographic units.

The protective clothing is particularly convenient to put on because it can be individually adjusted to the girth of the respective patient with a flick of the wrist. The front-mounted fastening system prevents parts of the protective clothing from entering the image-relevant area due to it being closed inaccurately. The back lies smooth and flat against the body.

- **Protective Material** Novalite Plus or Leadfree 2.0
- **Lead Equivalent** 0.50 mm Pb (Front and Back)
- **Size**

Adult	Front: 80.0 x 60.0 cm (H x W)
Children	Front: 60.0 x 45.0 cm (H x W)
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Hybrid, available in the colours:
Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy (ComforTex® HPMF only)

MAVIG GmbH

X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com