

EN

RD642

MAVIG



Models: **RD6425E** Adult patient X-ray apron for X-rays in a seated position
RD6425K Children patient X-ray apron for X-rays in a seated position

Made in Germany

CE 0302

RD642 – Patient X-Ray Apron for X-Rays in a Seated Position

The MAVIG patient protective apron RD642, which is used for dental X-rays in a seated position, is flexible and comfortable to wear.

The high-quality hook-and-loop fasteners are conveniently placed in the shoulder area, making it precise, secure and quick to put on and remove.

- **Protective Material** Novalite Plus or Leadfree 2.0
- **Lead Equivalent** 0.50 mm Pb (front)
- **Size**

| | |
|----------|-------------------------------|
| Adult | Front: 80.0 x 60.0 cm (H x W) |
| Children | Front: 60.0 x 45.0 cm (H x W) |
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Hybrid, available in the colours:
Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy (ComforTex® HPMF only)

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com