

EN

RD644

MAVIG



CAN BE USED ON BOTH SIDES

The model RD644 is designed as a reversible X-ray apron for patients.



Models: RD6445E
RD6445K

Adult patient X-ray apron with thyroid protection, for X-rays in a seated position
Children patient X-ray apron with thyroid protection, for X-rays in a seated position

Made in Germany

CE 0302

RD644 – Patient X-Ray Apron with Thyroid Protection

Just like the classic dental apron (RD642), this model offers high flexibility and wearing comfort. In addition to the large-area protection up to the gonads, the thyroid gland is also covered when taking X-rays of the seated patient.

The practical hook-and-loop fastener at the back of the neck on the thyroid protection collar ensures a snug, secure and quick fit. Loops on the back allow easy attachment to the dentist's chair.

- **Protective Material** Novalite Plus or Leadfree 2.0
- **Lead Equivalent** 0.50 mm Pb (front)
- **Size**

Adult	Front: 80.0 x 60.0 cm (H x W)
Children	Front: 60.0 x 45.0 cm (H x W)
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Hybrid, available in the colours:

<i>Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy (ComforTex® HPMF only)</i>	and two-colour combinations: <i>Curacao with Orchidee or Blue Admiral with Baccara</i>
--	---

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com