



Model: RP2705 Breast protection for female patients during mammographic radiation procedures

Made in Germany


CE 0302

RP270 – Breast Protection for Women

The mammary gland shield efficiently protects your patients during mammograms.

The shield consists of two equal parts, each 18 x 40 cm (W x H), connected by a lengthadjustable hook-and-loop fastener for comfortable donning and easy adjustment.

- **Protective Material** Standard Lead
- **Lead Equivalent** 0.50 mm Pb
- **Size** Universal
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Plus (design Happy only), available in the colours:
Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or designs: Happy
- **Edging** All edges are covered with a seem band.



MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com