

EN

RP274

MAVIG



Art.-No.: RP2741G6212 Gonad Protection Pants Set, 3 Children Sizes

Made in Germany

CE 0302

RP274 – Gonad Protection Pants Set for Children

The gonad protection pants comfortably ensure complete protection of the gonads of your smallest patients. The practical hook-and-loop fastening makes them quick and easy to put on and remove and enable convenient size adjustment.

The complete size set includes one pair of pants each for babies, toddlers and children from 1 to approx. 6 years of age.

- **Protective Material** Standard Lead
- **Lead Equivalent** 1.00 mm Pb
- **Outer Cover/Colour** ComforTex® HPMF Plus, in design Happy
- **Sizes** The 3-piece complete set RP2741G6212 is suitable for infants, toddlers and children from 1 to 6 years of age.

The respective sizes of each set are variably adjustable by means of hook-and-loop fastener.

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com