

EN

RP664

MAVIG



<b>Models:</b>	<b>RP66452S</b>	Children Protective Apron (SMALL)
	<b>RP66452M</b>	Children Protective Apron (MEDIUM)
	<b>RP66452L</b>	Children Protective Apron (LARGE)

Made in Germany

CE 0302

## Basic<sup>Mini</sup> – RP664 – Children's Front Protection Apron

This children's apron is especially designed for universal use and therefore the optimal front protection for your little patients. The generously cut, overlapping back panels provide a balanced fit.

Children's X-ray protection in particular, should be easy to put on and take off. The model RP664 can be easily closed in the front with a variable snap buckle closure leaving enough room for adjustments for an exact fit.

Children's X-ray protection is manufactured by MAVIG with the innovative Leadfree 2.0 radiation protection material to achieve the best possible weight ratio and thus impose only the most necessary physical strain on the children.

- **Protective Material** Leadfree 2.0
- **Lead Equivalent** Front protection 0.50 mm Pb / Back panels 0.25 mm Pb
- **Outer Cover/Colour** ComforTex® HPMF, available in the colours:  
*Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy*
- **Sizes**

SMALL	Length: 59 cm, Width: 35 cm, for children 3 to 5 years old
MEDIUM	Length: 67 cm, Width: 43 cm, for children 6 to 8 years old
LARGE	Length: 84 cm, Width: 45 cm, for children 9 to 12 years old

**MAVIG GmbH**  
X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone +49 (0) 89 / 420 96 - 0  
Fax +49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com